



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
10/11 Epoxide	Serum	Carbamazepine Epoxide
11 Deoxycortisol	Serum or Plasma (EDTA)	11-Deoxycortisol
11-Deoxycortisol	Serum or Plasma (EDTA)	
17 Ketosteroid	Not Available - See Note.	17-Ketosteroid
17 OH Corticosteroid	Not Available - See Note.	17-Hydroxycorticosteroid
17 Oh Progesterone	Plasma (Li-heparin)	17-Hydroxyprogesterone
17-Hydroxycorticosteroid	Not Available - See Note.	
17-Hydroxyprogesterone	Plasma (Li-heparin)	
17-KGS	Not Available - See Note.	17-Ketogenic Steroids
17-KS	Not Available - See Note.	17-Ketosteroid
17-Ketogenic Steroids	Not Available - See Note.	
17-Ketosteroid	Not Available - See Note.	
17-Oh Progesterone	Plasma (Li-heparin)	17-Hydroxyprogesterone
17-b-Estradiol	Serum or Plasma (Hep. or EDTA)	Estradiol
5' Nucleotidase	Not Available - See Note.	
5-Fluorocytosine	Not Available - See Note.	
5-HT	Serum	Serotonin
5-HT	Urine (24h)	Serotonin
5-Hydroxyindole Acetic Acid	Urine (24h)	
5-Hydroxyindole Acetic Acid	Urine (Random)	
5-Hydroxytryptamine	Serum	Serotonin
5-Hydroxytryptamine	Urine (24h)	Serotonin
5HIAA	Urine (24h)	5-Hydroxyindole Acetic Acid
5HIAA	Urine (Random)	5-Hydroxyindole Acetic Acid
5NT	Not Available - See Note.	5' Nucleotidase
6-Monoacetylmorphine	Urine (Random)	Heroin: Screen
7 Dehydrocholesterol	Amniotic Fluid	7-Dehydrocholesterol
7-DHC	Amniotic Fluid	7-Dehydrocholesterol
7-Dehydrocholesterol	Amniotic Fluid	
7-Dehydrocholesterol	Serum or Plasma (Heparin)	
7DHC	Serum or Plasma (Heparin)	7-Dehydrocholesterol
8-Hydroxyamoxapine	Not Available - See Note.	Amoxapine
8-Hydroxyloxapine	Not Available - See Note.	Loxapine
α-1-Acid Glycoprotein	Serum	
α-1-Antitrypsin	Serum	
α-1-Antitrypsin: Phenotype	Serum	
α-2-Antiplasmin	Plasma (Citrate)	
α-2-Macroglobulin	Serum	
α-Fetoprotein	Amniotic fluid	
α-Fetoprotein	Fluid	
α-Fetoprotein	Serum	
α-Glycoprotein Subunit	Serum	
α-Tocopherol	Serum	Vitamin E



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
β1C Globulin	Serum	C3 Complement
β1E Globulin	Serum	C4 Complement
β-2-Glycoprotein 1 Antibody	Serum	
β-2-Microglobulin	Serum	
β-2-Microglobulin	Urine (Fetal)	
β-2-Microglobulin	Urine (Random)	
β-2-Microglobulin - Fluid	Not available - See Note.	
β-Carotene	Serum	
β-Hexosaminidase	Serum	
β-Hydroxybutyrate	Serum or Plasma (Heparin or Fluoride)	
β-Transferrin	Fluid	
γ-Glutamyltransferase	Serum	
γ-Hydroxybutyrate	Urine (random)	
Δ-ALA	Urine (Random)	Porphyrin Precursors (ALA & PBG)
A2 Macroglobulin	Serum	α-2-Macroglobulin
α-Glycoprotein Subunits	Serum	α-Glycoprotein Subunit
ACE	CSF	Angiotensin Converting Enzyme
ACE	Serum	Angiotensin Converting Enzyme
ACP	Not Available - See Note.	Acid Phosphatase: Total
ACTH	Plasma (EDTA)	Adrenocorticotrophic Hormone
ADH	Plasma (EDTA)	Antidiuretic Hormone
ADIPOQ	Serum	Adiponectin
AFP	Amniotic fluid	α-Fetoprotein
AFP	Fluid	α-Fetoprotein
AFP	Serum	α-Fetoprotein
ALA	Urine (Random)	Porphyrin Precursors (ALA & PBG)
ALP	Serum	Alkaline Phosphatase
ALP	Serum	Alkaline Phosphatase: Isoenzymes
ALP	Not Available - See note.	Alkaline Phosphatase: Neutrophil
ALT	Serum	Alanine Aminotransferase
ANA	Fluid	Nuclear Antibodies
ANA	Serum	Nuclear Antibodies
ANCA	Serum	Neutrophil Cytoplasmic Antibody
ANF	Fluid	Nuclear Antibodies
ANF	Serum	Nuclear Antibodies
ANNA-1	CSF	Paraneoplastic Autoantibody Panel
ANNA-1	Serum	Paraneoplastic Autoantibody Panel
ANNA-2	CSF	Paraneoplastic Autoantibody Panel
ANNA-2	Serum	Paraneoplastic Autoantibody Panel
APCR	Plasma (Citrate)	Activated Protein C Resistance
ASA	Gastric lavage	Acetylsalicylic Acid
ASA	Serum	Acetylsalicylic Acid
ASA	Not Available - See Note.	Acetylsalicylic Acid (Urine)



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
ASA	Not Available - See Note.	Percodan
ASCA	Serum	Saccharomyces cer. Antibodies
ASOT	Not Available - See Note.	Anti-Streptolysin O Screen
ASOT	Serum	Anti-Streptolysin O Titre
AST	Serum	Aspartate Aminotransferase
Acetaminophen	Gastric lavage.	
Acetaminophen	Serum	
Acetaminophen	Not Available - See Note.	Percocet
Acetone	Gastric lavage	Alcohol Fract. and Quant.
Acetone	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Acetone	Urine (Random)	Alcohol Fract. and Quant.
Acetone	Vitreous Humour	Alcohol Fract. and Quant.
Acetylcholine Receptor Antibody	Serum	
Acetylcholinesterase: RBC	Not Available - See Note.	
Acetylsalicylic Acid	Gastric lavage	
Acetylsalicylic Acid	Serum	
Acetylsalicylic Acid	Not Available - See Note.	Percodan
Acetylsalicylic Acid (Urine)	Not Available - See Note.	
Acid Glycoprotein	Serum	α-1-Acid Glycoprotein
Acid Phosphatase: Prostatic	Not Available - See Note	
Acid Phosphatase: Total	Not Available - See Note.	
Activated Protein C Resistance	Plasma (Citrate)	
Acylcarnitine	Serum	
Adenovirus	Multiple	
Adenovirus Identification by PCR	Multiple	
Adiponectin	Serum	
Adrenal Antibodies	Serum	
Adrenaline	Plasma (EDTA)	Catecholamines
Adrenaline	Urine (24h)	Catecholamines
Adrenocorticotrophic Hormone	Plasma (EDTA)	
Ag	Erythrocyte	Silver
Ag	Urine (24h)	Silver
Ag	Urine (Random)	Silver
Al	Fluid	Aluminium
Al	Plasma (EDTA)	Aluminium
Al	Urine (24h)	Aluminium
Al	Urine (Random)	Aluminium
Alanine Aminotransferase	Serum	
Albumin	Fluid / Faeces	
Albumin	Serum	
Albumin	Urine (24h)	
Albumin	Urine (Random)	
Alcohol	Urine (Random)	Ethanol



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Alcohol Fract. and Quant.	Gastric lavage	
Alcohol Fract. and Quant.	Serum or Plasma (Fluoride)	
Alcohol Fract. and Quant.	Urine (Random)	
Alcohol Fract. and Quant.	Vitreous Humour	
Alcohols	Gastric lavage	Alcohol Fract. and Quant.
Alcohols	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Alcohols	Urine (Random)	Alcohol Fract. and Quant.
Alcohols	Vitreous Humour	Alcohol Fract. and Quant.
Aldactone	Urine (Random)	Spironolactone: Screen
Aldolase	Not Available - See Note.	
Aldosterone	Serum or Plasma (Hep)	
Aldosterone	Urine (24h)	
Alkaline Phosphatase	Serum	
Alkaline Phosphatase: Bone Specific	Plasma (Heparin) or Serum	
Alkaline Phosphatase: Isoenzymes	Serum	
Alkaline Phosphatase: Neutrophil	Not Available - See note.	
Allergens	Serum	Immunoglobulin E Ab. Single Allergen
Allergic Alveolitis (Lung Serology)	Serum	
Allergic Lung	Serum	Allergic Alveolitis (Lung Serology)
Alpha 1 Acid Glycoprotein	Serum	α -1-Acid Glycoprotein
Alpha 1 Antitrypsin	Serum	α -1-Antitrypsin
Alpha 1 Antitrypsin Phenotype	Serum	α -1-Antitrypsin: Phenotype
Alpha 2 Antiplasmin	Plasma (Citrate)	α -2-Antiplasmin
Alpha Glycoprotein Subunit	Serum	α -Glycoprotein Subunit
alpha Hydroxyprogesterone	Plasma (Li-heparin)	17-Hydroxyprogesterone
Alpha Protease Inhibitor	Serum	α -1-Antitrypsin: Phenotype
Alpha-2-Macroglobulin	Serum	α -2-Macroglobulin
Alpha-Fetoprotein	Amniotic fluid	α -Fetoprotein
Alpha-Fetoprotein	Fluid	α -Fetoprotein
Alpha-Fetoprotein	Serum	α -Fetoprotein
Alprazolam	Not Available - See Note.	
Aluminium	Fluid	
Aluminium	Plasma (EDTA)	
Aluminium	Urine (24h)	
Aluminium	Urine (Random)	
Aluminum	Fluid	Aluminium
Aluminum	Plasma (EDTA)	Aluminium
Aluminum	Urine (24h)	Aluminium
Aluminum	Urine (Random)	Aluminium
Amikacin: Peak	Serum	
Amikacin: Trough	Serum	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Amikin	Serum	Amikacin: Peak
Amikin	Serum	Amikacin: Trough
Amino Acids	CSF	Glycine
Amino Acids - Frac. & Quant. (Fluid)	Not Available - See Note.	
Amino Acids: Quantitation	Serum or Plasma (Hep. / EDTA)	
Amino Acids: Quantitation	Urine (Random)	
Aminoglycoside	Serum	Amikacin: Peak
Aminoglycoside	Serum	Amikacin: Trough
Aminoglycoside	Not Available - See Note.	Netilmicin
Aminoglycoside	Serum	Tobramycin: Trough
Aminoglycoside	Serum	Vancomycin: Peak
Aminoglycosides	Serum	Gentamicin: Peak
Aminoglycosides	Serum	Gentamicin: Trough
Aminoglycosides	Serum	Tobramycin: Peak
Aminoglycosides	Serum	Vancomycin: Trough
Aminolevulinic Acid	Urine (24h)	Porphyrin Precursors (ALA & PBG)
Aminolevulinic Acid	Urine (Random)	Porphyrin Precursors (ALA & PBG)
Amiodarone	Serum	
Amitriptyline	Gastric lavage	
Amitriptyline	Serum or plasma (EDTA)	
Amitriptyline	Urine (Random)	
Ammonia	Plasma (EDTA)	
Amobarbital	Gastric lavage	
Amobarbital	Serum	
Amoxapine	Not Available - See Note.	
Amphetamine	Urine (Random)	Methylenedioxy Amphetamine: Screen
Amphetamine	Urine (Random)	
Amylase	Serum and Urine	Amylase/Creatinine Clearance Ratio
Amylase	Fluid	
Amylase	Serum	
Amylase	Urine (Random or 24h)	
Amylase - Macro	Serum	Macroamylase
Amylase Isoenzyme	Serum	
Amylase/Creatinine Clearance Ratio	Serum and Urine	
Amyloidosis	Serum	Immunoglobulins: Free Light Chains
Amytal	Gastric lavage	Amobarbital
Amytal	Serum	Amobarbital
Anafranil	Gastric lavage	Clomipramine
Anafranil	Serum or plasma (EDTA)	Clomipramine
Anafranil	Urine (Random)	Clomipramine
Analogue-Tracer Free	Not Available - See Note.	Testosterone: Free

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Testosterone		
Ancotil	Not Available - See Note.	5-Fluorocytosine
Androstenedione	Serum	
Angel Dust	Urine (Random)	Phencyclidine: Screen
Angiotensin Converting Enzyme	CSF	
Angiotensin Converting Enzyme	Serum	
Anti GAD	Serum	Glutamic Acid Decarboxylase Antibodies
Anti-Adrenal	Serum	Adrenal Antibodies
Anti-CCP	Serum	Cyclic Citrullinated Peptide Antibodies
Anti-Cardiolipin	Serum	Cardiolipin Antibodies
Anti-DNA	Serum	DNA Antibodies
Anti-ENA	Serum	Extractable Nuclear Antigens Ab.
Anti-Endomysial	Serum	Endomysial Antibodies
Anti-GBM	Serum	Glomerular Base. Membrane Antibody
Anti-Glomerular Basement Membrane	Serum	Glomerular Base. Membrane Antibody
Anti-HAV	Multiple	Hepatitis A Antibody - IgG
Anti-HAV	Multiple	Hepatitis A Antibody - IgM
Anti-HBSAg	Multiple	Hepatitis B Surface Antigen Antibody
Anti-HBc	Multiple	Hepatitis B Core Antibody
Anti-HBe	Multiple	Hepatitis B e Antibody
Anti-HU	CSF	Paraneoplastic Autoantibody Panel
Anti-Heart Muscle	Serum	Heart Muscle Antibodies
Anti-Histones	Serum	Histone Antibodies
Anti-Hu	Serum	Paraneoplastic Autoantibody Panel
Anti-Insulin	Serum	Insulin Antibodies
Anti-Intrinsic Factor	Serum	Intrinsic Factor Antibodies
Anti-Islet Cell	Serum	Pancreatic Islet Cell Antibodies
Anti-LKM	Serum	Liver Kidney Microsomal Antibodies
Anti-Liver Kidney Microsomal	Serum	Liver Kidney Microsomal Antibodies
Anti-Microsomal	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Anti-Microsomal	Serum or Plasma (Na-Heparin)	Thyroid Peroxidase Antibodies
Anti-Milk	Not Available - See Note.	Milk Antibodies
Anti-Mitochondrial	Serum	Mitochondrial Antibodies
Anti-Mullerian Hormone	Serum or Plasma (EDTA or Hep)	
Anti-Myocardial	Serum	Heart Muscle Antibodies
Anti-Nuclear Factor	Serum	Nuclear Antibodies
Anti-Ovarian	Not Available - See Note.	Ovarian Antibodies
Anti-PCNA	Serum	Proliferating Cell Nuclear Antibodies
Anti-Pancreatic Islet Cell	Serum	Pancreatic Islet Cell Antibodies

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Anti-Parietal	Serum	Parietal Antibodies
Anti-Parotid	Not Available - See Note.	Parotid Antibodies
Anti-Peroxidase	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Anti-Platelet	Not Available - See Note.	Platelet Associated IgG
Anti-RI	CSF	Paraneoplastic Autoantibody Panel
Anti-Ri	Serum	Paraneoplastic Autoantibody Panel
Anti-Skeletal	Serum	Skeletal Muscle Antibodies
Anti-Skin	Serum	Skin Antibodies
Anti-Smooth Muscle	Serum	Smooth Muscle Antibodies
Anti-Streptolysin O Screen	Not Available - See Note.	
Anti-Streptolysin O Titre	Serum	
Anti-Striated Muscle	Serum	Skeletal Muscle Antibodies
Anti-Thrombin III: Activity	Plasma (Citrate)	
Anti-Thrombin III: Antigen	Plasma (Citrate)	
Anti-Thrombin III: Functional	Plasma (Citrate)	Anti-Thrombin III: Activity
Anti-Thyroglobulin	Serum or Plasma (Na-Heparin)	Thyroglobulin Antibodies
Anti-Thyroglobulin	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Anti-Thyroid	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Anti-Thyrotropin Receptor	Serum	Thyroid Receptor Antibody
Anti-YO	CSF	Paraneoplastic Autoantibody Panel
Anti-Yo	Serum	Paraneoplastic Autoantibody Panel
Anti-freeze	Serum or Plasma (Fluoride)	Ethylene Glycol
Anti-freeze	Urine (Random)	Ethylene Glycol
Anti-freeze	Gastric lavage	Ethylene Glycol
Anticonvulsant	Serum	Lamotrigine
Antidepressant	Urine (Random)	Amphetamine
Antidepressant	Gastric lavage	Carbamazepine
Antidepressant	Serum	Carbamazepine
Antidepressant	Serum	Fluoxetine
Antidepressant	Serum	Fluvoxamine
Antidepressant	Serum	Lamotrigine
Antidepressant	Serum	Lithium
Antidepressant	Urine (Random)	Methylphenidate: Screen
Antidepressant	Serum or Plasma (Heparin)	Paroxetine
Antidepressant	Gastric lavage	Protriptyline
Antidepressant	Serum or plasma (EDTA)	Protriptyline
Antidepressant	Urine (Random)	Protriptyline
Antidepressant	Serum or plasma (Heparin)	Sertraline
Antidepressant	Serum	Trazodone
Antidepressant	Serum	Valproic Acid
Antidepressant	Serum or Plasma (Heparin)	Valproic Acid: Free



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Antidepressant	Urine (Random)	Venlafaxine: Screen
Antidiuretic Hormone	Plasma (EDTA)	
Antihaemophilic Factor	Plasma (Citrate)	Coagulation Factor VIII
Antihaemophilic Factor	Plasma (Citrate)	Von Willebrand Factor Antigen
Antihaemophilic Factor	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Antimony	Urine (24h)	
Antimony	Urine (Random)	
Antimony	Whole blood (Heparin)	
Antimony - Plasma	Not Available - See Note.	
Antinuclear Antibody	Fluid	Nuclear Antibodies
Antinuclear Antibody	Serum	Nuclear Antibodies
Antiplasmin	Plasma (Citrate)	α -2-Antiplasmin
Antithrombin	Plasma (Citrate)	Anti-Thrombin III: Activity
Antithrombin	Plasma (Citrate)	Anti-Thrombin III: Antigen
Antithrombin III: Immunological	Plasma (Citrate)	Anti-Thrombin III: Antigen
Antitrypsin	Serum	α -1-Antitrypsin
Antitrypsin	Serum	α -1-Antitrypsin: Phenotype
Apo B	Serum or Plasma (EDTA)	Apolipoprotein B
Apo C	Plasma (EDTA)	Apolipoprotein C &/or E
Apo E	Plasma (EDTA)	Apolipoprotein C &/or E
Apolipoprotein A1	Serum or Plasma (EDTA)	
Apolipoprotein B	Serum or Plasma (EDTA)	
Apolipoprotein C &/or E	Plasma (EDTA)	
Apolipoprotein C2 Activation	Plasma (EDTA)	
Aquaporin 4 Antibody	Serum	Neuromyelitis Optica Antibody (IgG)
Arginine	Urine (Random)	Cystinuria Monitoring
Arginine Vasopressin	Plasma (EDTA)	Antidiuretic Hormone
Arixtra	Plasma (Citrate)	Heparin: Fondaparinux
Arsenic	Whole Blood (EDTA)	
Arsenic	Hair	
Arsenic	Nail	
Arsenic: Inorganic	Urine (24h)	
Arsenic: Inorganic	Urine (Random)	
Arsenic: Total	Urine (24h)	
Arsenic: Total	Urine (Random)	
Aryl Sulfatase A	Whole blood (Heparin)	
As	Whole Blood (EDTA)	Arsenic
As	Hair	Arsenic
As	Nail	Arsenic
As	Urine (24h)	Arsenic: Inorganic
As	Urine (Random)	Arsenic: Inorganic
As	Urine (24h)	Arsenic: Total
As	Urine (Random)	Arsenic: Total
Ascendin	Not Available - See Note.	Amoxapine



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Ascorbic Acid	Plasma (Heparin)	Vitamin C
Aspartate Aminotransferase	Serum	
Aspergillus	Serum	Allergic Alveolitis (Lung Serology)
Aspirate	Aspirate	Aspiration Biopsy
Aspirates	Variable	Pathology Services
Aspiration Biopsy	Aspirate	
Aspirin	Gastric lavage	Acetylsalicylic Acid
Aspirin	Serum	Acetylsalicylic Acid
Ativan	Gastric lavage	Lorazepam
Ativan	Serum	Lorazepam
Ativan	Urine (Random)	Lorazepam
Au	Serum or Plasma (EDTA)	Gold
Aureobasidium	Serum	Allergic Alveolitis (Lung Serology)
Australia Antigen	Multiple	Hepatitis B Surface Antigen
Aventyl	Gastric lavage	Nortriptyline
Aventyl	Serum or Plasma (EDTA)	Nortriptyline
Aventyl	Urine (Random)	Nortriptyline
Avian Precipitins: Budgie	Serum	
Avian Precipitins: Canary	Serum	
Avian Precipitins: Chicken	Serum	
Avian Precipitins: Cockatiel	Serum	
Avian Precipitins: Duck	Serum	
Avian Precipitins: Goose	Serum	
Avian Precipitins: Parrot	Serum	
Avian Precipitins: Turkey	Serum	
Avian precipitins: Pigeon	Serum	
B-CTx	Plasma (EDTA)	C-Telopeptide
B-Cells	Whole blood (EDTA)	Lymphocyte Marker Analysis
B-Type Natriuretic Peptide (NT-Pro)	Serum	
B2MG	Serum	β-2-Microglobulin
B2MG	Urine (Fetal)	β-2-Microglobulin
B2MG	Urine (Random)	β-2-Microglobulin
BK Virus	Plasma (EDTA)	
BK Virus	Urine (Random)	
BUN	Serum	Urea
BUN	Urine (24h)	Urea
Barbital	Gastric lavage	
Barbital	Serum	
Barbitone	Serum	Barbital
Barbiturates: Screen	Serum	
Barbiturates: Screen	Urine (Random)	
Barium	Urine (24h)	
Barium	Urine (Random)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Benadryl	Serum	Diphenhydramine
Benzodiazepine	Serum	Chlordiazepoxide
Benzodiazepine	Gastric lavage	Clobazam
Benzodiazepine	Urine (Random)	Clobazam
Benzodiazepine	Serum	Clobazam
Benzodiazepine	Gastric lavage	Clonazepam
Benzodiazepine	Serum	Clonazepam
Benzodiazepine	Urine (Random)	Clonazepam
Benzodiazepine	Serum	Clorazepate
Benzodiazepine	Serum	Diazepam
Benzodiazepine	Urine (Random)	Flunitrazepam
Benzodiazepine	Serum	Lorazepam
Benzodiazepine	Urine (Random)	Lorazepam
Benzodiazepine	Gastric lavage	Nitrazepam
Benzodiazepine	Serum	Nitrazepam
Benzodiazepine	Urine (Random)	Nitrazepam
Benzodiazepine	Serum	Oxazepam
Benzodiazepine Identification	Urine (random)	
Benzodiazepines: Screen	Serum	
Benzodiazepines: Screen	Urine (Random)	
Benzoyllecgonine	Urine (Random)	Cocaine: Screen
Benzotropine Mesylate: Screen	Urine (Random)	
Beryllium	Urine (24h)	
Beryllium	Urine (Random)	
Beta 1E Globulin	Serum	C4 Complement
Beta 1c Globulin	Serum	C3 Complement
Beta 2 Microglobulin	Serum	β -2-Microglobulin
Beta 2 Microglobulin	Urine (Random)	β -2-Microglobulin
Beta Carotene	Serum	β -Carotene
Beta Hexosaminidase	Serum	β -Hexosaminidase
Beta Transferrin	Fluid	β -Transferrin
Beta-Crosslaps	Plasma (EDTA)	C-Telopeptide
Beta-Hydroxybutyrate	Serum or Plasma (Heparin or Fluoride)	β -Hydroxybutyrate
Bile Acids: Total	Serum	
Bile Salts	Serum	Bile Acids: Total
Bilirubin	Amniotic fluid	
Bilirubin: Direct	Serum	
Bilirubin: Total	Serum	
Bioavailable Testosterone	Serum	Testosterone: Bioavailable
Biotinidase	Serum or Plasma	
Biphenin	Urine (Random)	Methylphenidate: Screen
Biquin	Serum	Quinidine
Bird Fanciers	Serum	Allergic Alveolitis (Lung Serology)

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Bismuth	Urine (Random)	
Blastomycosis	Multiple	
Blood Urea Nitrogen	Serum	Urea
Blood Urea Nitrogen	Urine (24h)	Urea
Bone ALP	Plasma (Heparin) or Serum	Alkaline Phosphatase: Bone Specific
Boron	Urine (24h)	
Boron	Urine (Random)	
Bromazepam	Not Available - See Note.	
Bromide	Not Available - See Note.	
Brucellosis	Multiple	
Budgie Antigen	Serum	Avian Precipitins: Budgie
Bupivacaine	Not Available - See Note.	
Butabarbital	Gastric lavage	
Butabarbital	Serum	
Butabarbital	Urine	
Butalbital	Serum	
Butisol	Gastric lavage	Butabarbital
Butisol	Serum	Butabarbital
Butisol	Urine	Butabarbital
C-ANCA	Serum	Neutrophil Cytoplasmic Antibody (C only)
C-Peptide	Serum	
C-Reactive Protein	Serum.	
C-Telopeptide	Plasma (EDTA)	
C1 Esterase Inhibitor: Functional	Plasma (Citrate)	
C1 Esterase Inhibitor: Immunological	Serum	
C1 Function	Serum	
C1Q Binding	Serum	Immune Complexes
C1q Complement Component	Plasma (EDTA)	
C1q Level	Plasma (EDTA)	C1q Complement Component
C1r Complement Component	Plasma (EDTA)	
C1r Level	Plasma (EDTA)	C1r Complement Component
C1s	Plasma (EDTA)	C1s Complement Component
C1s Complement Component	Plasma (EDTA)	
C2 Complement	Serum	
C2 Function	Serum	
C2 Protein	Serum	C2 Complement
C3 Complement	Serum	
C4 Complement	Serum	
C4 Complement: Fluid	Not Available - See Note.	
C5 Complement	Plasma (EDTA)	
C5 Level	Plasma (EDTA)	C5 Complement
C6 Complement	Plasma (EDTA)	
C7 Complement	Plasma (EDTA)	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
C7 Level	Plasma (EDTA)	C7 Complement
C8 Complement	Plasma (EDTA)	
C8 Level	Plasma (EDTA)	C8 Complement
C9 Complement	Plasma (EDTA)	
C9 Level	Plasma (EDTA)	C9 Complement
CA 125	Fluid	
CA 125	Serum	
CA 15-3	Serum or Plasma (EDTA)	
CA 19-9	Fluid	
CA 19-9	Serum	
CA125	Fluid	CA 125
CA125	Serum	CA 125
CA15-3	Serum or Plasma (EDTA)	CA 15-3
CA153	Serum or Plasma (EDTA)	CA 15-3
CA19-9	Fluid	CA 19-9
CA19-9	Serum	CA 19-9
CA199	Fluid	CA 19-9
CA199	Serum	CA 19-9
CANCA	Serum	Neutrophil Cytoplasmic Antibody
CANCA	Serum	Neutrophil Cytoplasmic Antibody (C only)
CCP Antibody	Serum	Cyclic Citrullinated Peptide Antibodies
CD4	Whole blood (EDTA)	Lymphocyte Marker Analysis
CD8	Whole blood (EDTA)	Lymphocyte Marker Analysis
CEA	Fluid	Carcinoembryonic Antigen
CEA	Serum	Carcinoembryonic Antigen
CH50	Fluid (Synovial) and Serum	Complement: Total
CH50	Serum	Complement: Total
CIC	Serum	Immune Complexes
CK	Serum or Plasma (Hep. or EDTA)	Creatine Kinase
CK Isoenzyme	Serum	Creatine Kinase: CK-MB Isoenzyme
CK Total	Serum or Plasma (Hep. or EDTA)	Creatine Kinase
CK-2	Serum	Creatine Kinase: CK-MB Isoenzyme
CK2	Serum	Creatine Kinase: CK-MB Isoenzyme
CKMB	Serum	Creatine Kinase: CK-MB Isoenzyme
CKT	Serum or Plasma (Hep. or EDTA)	Creatine Kinase
CMV	Multiple	Cytomegalovirus
CO	Whole blood (Heparin)	Carbon Monoxide
CPK	Serum or Plasma (Hep. or EDTA)	Creatine Kinase
CRP	Serum.	C-Reactive Protein
Ca	Erythrocyte	Calcium



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Ca	Serum	Calcium
Ca	Urine (24 h)	Calcium
Ca	Serum	Calcium: Ionized
Cadmium	Erythrocyte	
Cadmium	Urine (24h)	
Cadmium	Urine (Random)	
Cadmium	Whole blood (EDTA)	
Caffeine	Serum	
Caffeine - Urine	Not Available - See Note.	
Caffeine Quantitation	Serum	Caffeine
Caffeine Screen	Not Available - See Note.	Caffeine - Urine
Calcidiol	Serum or plasma (Heparin)	Vitamin D (25-Hydroxy)
Calcitonin	Serum	
Calcitriol	Serum or Plasma (Heparin)	Vitamin D (1,25-Dihydroxy)
Calcium	Erythrocyte	
Calcium	Fluid	
Calcium	Serum	
Calcium	Urine (24 h)	
Calcium	Urine (Random)	
Calcium Pyrophosphate	Synovial fluid	Synovial Fluid Examination
Calcium: Ionized	Serum	
Calculus Analysis	Calculus	
California Encephalitis Virus	Multiple	
cAMP	Not Available - See Note.	Cyclic AMP: Total
Canary Antigen	Serum	Avian Precipitins: Canary
Cannabinoids: Screen	Urine (Random)	
Cannabis	Urine (Random)	Cannabinoids: Screen
Carbamazepine	Gastric lavage	
Carbamazepine	Serum	
Carbamazepine Epoxide	Serum	
Carbon Monoxide	Whole blood (Heparin)	
Carboxyhemoglobin	Whole blood (Heparin)	Carbon Monoxide
Carcinoembryonic Antigen	Fluid	
Carcinoembryonic Antigen	Serum	
Cardiolipin Antibodies	Serum	
Carnitine	Serum	
Carnitine	Urine (Random)	
Carotene	Serum	β-Carotene
Catecholamine Metabolite	Plasma (EDTA)	Metanephrines
Catecholamine Metabolite	Urine (24h)	Metanephrines
Catecholamine Metabolite	Urine (24h)	Vanillylmandelic Acid
Catecholamines	Urine (24h)	Dopamine
Catecholamines	Plasma (EDTA)	
Catecholamines	Urine (24h)	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Cd	Erythrocyte	Cadmium
Cd	Urine (24h)	Cadmium
Cd	Urine (Random)	Cadmium
Cd	Whole blood (EDTA)	Cadmium
Celiac	Serum	Celiac Profile
Celiac	Serum	Gliadin (Deamidated) Antibodies
Celiac	Serum	Reticulin Antibodies
Celiac	Serum	Tissue Transglutaminase Antibodies
Celiac Profile	Serum	
Cellcept	Plasma (EDTA)	Mycophenolic Acid
Celontin	Serum or Plasma (Heparin)	Methsuximide
Centromere Antibody	Serum	Nuclear Antibodies
Ceruloplasmin	Serum	
Cervical specimens	Variable	Pathology Services
Chicken Antigen	Serum	Avian Precipitins: Chicken
Chlamydia	Multiple	
Chlamydia Isolation	Multiple	
Chloral Hydrate	Not Available - See Note.	
Chloramphenicol	Serum	
Chlordiazepoxide	Serum	
Chloride	Fluid / Faeces	
Chloride	Serum or Plasma (Heparin)	
Chloride	Urine (24h)	
Chloromycetin	Serum	Chloramphenicol
Chloroquine	Not Available - See Note	
Chlorothiazide	Not Available - See Note.	
Chlorpromazine	Gastric lavage	
Chlorpromazine	Serum	
Chlorpromazine	Urine (Random)	
Chlorpropamide	Not Available - See Note.	
Cholesterol	Serum or Plasma (Hep. or EDTA)	
Cholinesterase: Phenotype	Serum	
Cholinesterase: Total Activity	Serum	
Chorionic Gonadotrophin	Serum or Plasma (Hep. or EDTA)	
Chorionic Gonadotrophin: Tumour Marker	Serum	
Christmas Factor	Plasma (Citrate)	Coagulation Factor IX
Chromium	Erythrocyte	
Chromium	Plasma (EDTA)	
Chromium	Urine (24h)	
Chromium	Urine (Random)	
Chromium	Whole blood (EDTA)	
Chromium	Fluid	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Chromogranin A	Plasma (EDTA)	
Chromosome Studies	Various	
Cimetidine: Screen	Urine (Random)	
Circulating Anticoagulant	Plasma (Citrate)	Non-Spec. Coagulation Inhib.
Circulating Immune Complexes	Serum	Immune Complexes
Citrate	Urine (24h)	Citric Acid
Citrate	Urine (Random)	Citric Acid
Citric Acid	Urine (24h)	
Citric Acid	Urine (Random)	
Cl	Fluid / Faeces	Chloride
Cl	Serum or Plasma (Heparin)	Chloride
Cl	Urine (24h)	Chloride
Clauss Fibrinogen	Plasma (Sodium Citrate)	Fibrinogen: Immunological
Clobazam	Urine (random)	Benzodiazepine Identification
Clobazam	Gastric lavage	
Clobazam	Urine (Random)	
Clobazam	Serum	
Clomipramine	Gastric lavage	
Clomipramine	Serum or plasma (EDTA)	
Clomipramine	Urine (Random)	
Clonazepam	Urine (random)	Benzodiazepine Identification
Clonazepam	Gastric lavage	
Clonazepam	Serum	
Clonazepam	Urine (Random)	
Clorazepate	Serum	Diazepam
Clorazepate	Serum	
Clostridium difficile	Faeces	
Clozapine	Plasma (EDTA or Heparin)	
Clozaril	Plasma (EDTA or Heparin)	Clozapine
Cn	Whole blood (EDTA)	Cyanide
Co	Plasma (EDTA)	Cobalt
Co	Urine (24h)	Cobalt
Co	Urine (Random)	Cobalt
Coag. Factor VIII Inhibitor: Porcine	Not Available - See Note.	
Coagulation Factor	Plasma (Citrate)	Von Willebrand Assessment
Coagulation Factor 2	Plasma (Citrate)	Coagulation Factor II: Ecarin
Coagulation Factor 2: Biological	Plasma (Citrate)	Coagulation Factor II: Biological
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor II: Biological
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor II: Ecarin
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor IX
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor V
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor VII
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor VIII



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor X
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor XI
Coagulation Factor Assay	Plasma (Citrate)	Coagulation Factor XII
Coagulation Factor Assay	Plasma (Citrate)	Von Willebrand Factor Antigen
Coagulation Factor Assay	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Coagulation Factor II: Biological	Plasma (Citrate)	
Coagulation Factor II: Ecarin	Plasma (Citrate)	
Coagulation Factor IX	Plasma (Citrate)	
Coagulation Factor V	Whole Blood (ACD)	Factor V Leiden Mutation
Coagulation Factor V	Plasma (Citrate)	
Coagulation Factor VII	Plasma (Citrate)	
Coagulation Factor VIII	Plasma (Citrate)	
Coagulation Factor VIII Inhibitor: Human	Plasma (Citrate)	
Coagulation Factor VIII: Biological	Plasma (Citrate)	Coagulation Factor VIII
Coagulation Factor X	Plasma (Citrate)	
Coagulation Factor XI	Plasma (Citrate)	
Coagulation Factor XII	Plasma (Citrate)	
Coagulation Factor XIII: Screen	Plasma (Citrate)	
Cobalt	Plasma (EDTA)	
Cobalt	Plasma (K2 EDTA)	
Cobalt	Urine (24h)	
Cobalt	Urine (Random)	
Cocaine: Screen	Urine (Random)	
Coccidiomycosis	Multiple	
Cockatiel Antigen	Serum	Avian Precipitins: Cockatiel
Codeine	Urine (Random)	Narcotic Analgesics
Codeine: Screen	Urine (Random)	
Coeliac	Serum	Celiac Profile
Coeliac	Serum	Gliadin (Deamidated) Antibodies
Coeliac	Serum	Tissue Transglutaminase Antibodies
Cogentin	Urine (Random)	Benzotropine Mesylate: Screen
Cold Agglutinins	Serum	
Cold Hemolysis	Serum	Cold Agglutinins
Complement: Total	Fluid (Synovial) and Serum	
Complement: Total	Serum	
Compound S	Serum or Plasma (EDTA)	11-Deoxycortisol
Conjugated Bilirubin	Serum	Bilirubin: Direct
Copper	Erythrocyte	
Copper	Fluid	
Copper	Plasma (EDTA)	
Copper	Urine (24h)	
Copper	Urine (Random)	
Copper	Whole Blood (K2 EDTA)	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Copper	Tissue	
Coproporphyrin	Faeces (Random)	Porphyryns: Quantitation
Coproporphyrin	Urine (24h)	Porphyryns: Quantitation
Coproporphyrin	Urine (Random)	Porphyryns: Quantitation
Cordarone	Serum	Amiodarone
Cortisol	Saliva	
Cortisol	Serum	
Cortisol	Urine (24h)	
Cotinine	Urine (Random)	Nicotine: Screen
Coumadin	Plasma (Heparin)	Warfarin
Coxsackie A9	Multiple	
Coxsackie B1	Multiple	
Coxsackie B2	Multiple	
Coxsackie B3	Multiple	
Coxsackie B4	Multiple	
Coxsackie B5	Multiple	
Coxsackie B6	Multiple	
Cr	Erythrocyte	Chromium
Cr	Plasma (EDTA)	Chromium
Cr	Urine (24h)	Chromium
Cr	Urine (Random)	Chromium
Cr	Whole blood (EDTA)	Chromium
Cr	Fluid	Chromium
Creatine Kinase	Serum or Plasma (Hep. or EDTA)	
Creatine Kinase: CK-MB Isoenzyme	Serum	
Creatinine	Fluid	
Creatinine	Serum or Plasma (Heparin)	
Creatinine	Urine (24h)	
Creatinine (Amniotic Fluid)	Not Available - See Note.	
Creatinine Clearance	Serum and Urine (24h)	
Cryofibrinogen	Serum & Plasma (EDTA)	
Cryoglobulins	Serum	
Cryptococcosis	Multiple	
Crystal Identification	Synovial fluid	Synovial Fluid Examination
Cu	Erythrocyte	Copper
Cu	Fluid	Copper
Cu	Plasma (EDTA)	Copper
Cu	Urine (24h)	Copper
Cu	Urine (Random)	Copper
Cu	Whole Blood (K2 EDTA)	Copper
Cu	Tissue	Copper
Cyanide	Whole blood (EDTA)	
Cyanocobalamin	Serum	Vitamin B12



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Cyclic AMP: Total	Not Available - See Note.	
Cyclic Citrullinated Peptide Antibodies	Serum	
Cyclobenzaprine	Urine (Random)	
Cyclosporine	Whole Blood (EDTA)	
Cystatin C	Serum	
Cystine	Urine (Random)	Cystinuria Monitoring
Cystinuria Monitoring	Urine (Random)	
Cytology	Aspirate	Aspiration Biopsy
Cytology	Fluid	
Cytology	Sputum	
Cytomegalovirus	Multiple	
D-Dimer	Plasma (Citrate)	
D-Lactate	Serum or plasma (EDTA)	D-Lactic Acid
D-Lactic Acid	Serum or plasma (EDTA)	
DHEAS	Serum	Dehydroepiandrosterone Sulfate
DHT	Serum	Dihydrotestosterone
DNA Analysis	Various	Molecular Genetics
DNA Antibodies	Serum	
Dalmane	Not Available - See Note.	Flurazepam
danaparoid sodium	Plasma (Citrate)	Heparin: Orgaran Heparinoid
Darvon	Gastric lavage	Propoxyphene
Darvon	Urine (Random)	Propoxyphene
Darvon	Serum	Propoxyphene
Deamidated Gliadin Antibodies	Serum	Celiac Profile
Dehydroepiandrosterone Sulfate	Serum	
Demerol	Urine (Random)	Meperidine: Screen
Dengue	Multiple	
Deoxyypyridinoline	Not available - See Note.	
Depakene	Serum	Valproic Acid
Depakene	Serum or Plasma (Heparin)	Valproic Acid: Free
Desalkylflurazepam	Not Available - See Note.	Flurazepam
Desethylamiodarone	Serum	Amiodarone
Desipramine	Gastric lavage	Imipramine
Desipramine	Serum or plasma (EDTA)	Imipramine
Desipramine	Urine (Random)	Imipramine
Desipramine	Gastric lavage	
Desipramine	Serum or plasma (EDTA)	
Desipramine	Urine (Random)	
Desmethylclomipramine	Serum or plasma (EDTA)	Clomipramine
Desmethylclomipramine	Urine (Random)	Clomipramine
Desmethyldiazepam	Serum	Diazepam
Desyrel	Serum	Trazodone
Deuteroporphyrin	Faeces (Random)	Porphyryns: Quantitation

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Dextromethorphan	Urine or gastric lavage	
Diastase	Serum	Amylase
Diazepam	Serum	
Differential Count	Synovial fluid	Synovial Fluid Examination
Digoxin	Serum	
Digoxin - Free	Serum	Digoxin: Free
Digoxin: Free	Serum	
Dihydrotestosterone	Serum	
Dilantin	Serum	Phenytoin
Dilantin	Serum	Phenytoin: Free
Dimenhydrinate	Serum	
Diphenhydramine	Serum	
Diphenylhydantoin	Serum	Phenytoin: Free
Diphtheria Anti-toxin	Multiple	
Disopyramide	Plasma (Heparin)	
Diuril	Not Available - See Note.	Chlorothiazide
Dolophine	Fluid	Methadone
Dolophine	Serum	Methadone
Dolophine	Urine (Random)	Methadone
Dopamine	Urine (24h)	
Doriden	Gastric lavage	Glutethimide
Doriden	Serum	Glutethimide
Doriden	Urine (Random)	Glutethimide
Doxepin	Gastric lavage	
Doxepin	Serum or plasma (EDTA)	
Doxepin	Urine (Random)	
Drug Abuse Screen	Urine (Random)	
Drug Screen	Urine (Random)	Drug Abuse Screen
Drug Screen	Gastric lavage	
Drug Screen	Serum	
Drug Screen	Urine (Random)	
dsDNA	Serum	DNA Antibodies
Duck Antigen	Serum	Avian Precipitins: Duck
E1	Serum or Plasma (Heparin)	Estrone
E2	Serum or Plasma (Hep. or EDTA)	Estradiol
E3	Serum	Estriol: Unconjugated
EBV	Multiple	Epstein-Barr Virus
ECHIS	Plasma (Citrate)	Coagulation Factor II: Ecarin
EFT	Plasma (Citrate)	Euglobulin Fibrinolysis Time
EPO	Serum or Plasma (Heparin)	Erythropoietin
Eastern Equine Encephalitis Virus	Multiple	
Echinococcosis	Multiple	
Echovirus 1	Multiple	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Echovirus 2	Multiple	
Ecstasy	Urine (Random)	Methylenedioxy Amphetamine: Screen
Effexor	Urine (Random)	Venlafaxine: Screen
eGFR	Serum or Plasma (Heparin)	Creatinine
Elastase 1	Faeces	
Elavil	Gastric lavage	Amitriptyline
Elavil	Serum or plasma (EDTA)	Amitriptyline
Elavil	Urine (Random)	Amitriptyline
Eltroxin	Not Available - See Note.	Thyroxine (T4)
Endomysial Antibodies	Serum	
Enolase, Neuron Specific	Not Available - See Note.	
Entamoeba Histolytica	Multiple	
Epidermal Antibodies	Serum	Skin Antibodies
Epinephrine	Plasma (EDTA)	Catecholamines
Epinephrine	Urine (24h)	Catecholamines
Epival	Serum	Valproic Acid
Epival	Serum or Plasma (Heparin)	Valproic Acid: Free
Epstein-Barr Virus	Multiple	
Equanil	Gastric lavage	Meprobamate
Equanil	Serum	Meprobamate
Equanil	Urine (Random)	Meprobamate
Erythropoietin	Serum or Plasma (Heparin)	
Estimated Glomerular Filtration Rate	Serum or Plasma (Heparin)	Creatinine
Estradiol	Serum or Plasma (Hep. or EDTA)	
Estriol: Unconjugated	Serum	
Estrone	Serum or Plasma (Heparin)	
EtOH	Urine (Random)	Ethanol
Ethanol	Gastric lavage	Alcohol Fract. and Quant.
Ethanol	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Ethanol	Urine (Random)	Alcohol Fract. and Quant.
Ethanol	Vitreous Humour	Alcohol Fract. and Quant.
Ethanol	Urine (Random)	
Ethchlorvynol	Not Available - See Note.	
Ethosuximide	Serum	
Ethyl Alcohol	Urine (Random)	Ethanol
Ethyl Glucuronide	Urine (Random)	
Ethylene Glycol	Serum or Plasma (Fluoride)	
Ethylene Glycol	Urine (Random)	
Ethylene Glycol	Gastric lavage	
Euglobulin Clot Lysis	Plasma (Citrate)	Euglobulin Fibrinolysis Time
Euglobulin Fibrinolysis Time	Plasma (Citrate)	
Extacy	Urine (Random)	Methylenedioxy Amphetamine:

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
		Screen
Extasy	Urine (Random)	Methylenedioxy Amphetamine: Screen
Extractable Nuclear Antigens Ab.	Serum	
FDP	Plasma (Citrate)	D-Dimer
FFA	Serum	Fatty Acids, Free: Total
FK506	Whole Blood (EDTA)	Tacrolimus (FK506)
FSH	Serum or Plasma (Hep. or EDTA)	Follicle Stimulating Hormone
FT3	Serum or Plasma (Li-Hep./EDTA)	Triiodothyronine: Free
Factor 10	Plasma (Citrate)	Coagulation Factor X
Factor 11	Plasma (Citrate)	Coagulation Factor XI
Factor 12	Plasma (Citrate)	Coagulation Factor XII
Factor 13	Plasma (Citrate)	Coagulation Factor XIII: Screen
Factor 2	Plasma (Citrate)	Coagulation Factor II: Biological
Factor 2	Plasma (Citrate)	Coagulation Factor II: Ecarin
Factor 5	Plasma (Citrate)	Coagulation Factor V
Factor 7	Plasma (Citrate)	Coagulation Factor VII
Factor 8 Activity	Plasma (Citrate)	Coagulation Factor VIII
Factor 8 Antigen	Plasma (Citrate)	Von Willebrand Factor Antigen
Factor 8 Antigen	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Factor 8 C	Plasma (Citrate)	Coagulation Factor VIII
Factor 8 Functional	Plasma (Citrate)	Coagulation Factor VIII
Factor 8 Inhibitor	Plasma (Citrate)	Coagulation Factor VIII Inhibitor: Human
Factor 8 Related Antigen	Plasma (Citrate)	Von Willebrand Assessment
Factor 8 Related Antigen	Plasma (Citrate)	Von Willebrand Factor Antigen
Factor 8 Related Antigen	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Factor 9	Plasma (Citrate)	Coagulation Factor IX
Factor A	Serum	C3 Complement
Factor II	Plasma (Citrate)	Coagulation Factor II: Biological
Factor IX	Plasma (Citrate)	Coagulation Factor IX
Factor V	Plasma (Citrate)	Coagulation Factor V
Factor V Leiden Mutation	Whole Blood (ACD)	
Factor VII	Plasma (Citrate)	Coagulation Factor VII
Factor VIII	Plasma (Citrate)	Coagulation Factor VIII
Factor VIII	Plasma (Citrate)	Von Willebrand Assessment
Factor VIII	Plasma (Citrate)	Von Willebrand Factor Antigen
Factor VIII	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Factor X	Plasma (Citrate)	Coagulation Factor X
Factor XI	Plasma (Citrate)	Coagulation Factor XI
Factor XII	Plasma (Citrate)	Coagulation Factor XII
Farmers Lung Serology	Serum	Allergic Alveolitis (Lung Serology)



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Fat Globules	Faeces (Random)	Fat: Screen
Fat Screen	Faeces (Random)	Fat: Screen
Fat Screen	Urine (Random)	Fat: Screen
Fat: Quantitation	Faeces (72h)	
Fat: Screen	Faeces (Random)	
Fat: Screen	Urine (Random)	
Fatty Acids, Free: Total	Serum	
Fatty Acids, Long Chain	Serum or plasma (EDTA)	
Fe	Serum	Iron
Fe	Urine (24h)	Iron
Fe	Tissue	Iron
Fecal Elastase	Faeces	Elastase 1
Fecal Fat	Faeces (72h)	Fat: Quantitation
Fecal Fat Screen	Faeces (Random)	Fat: Screen
Ferritin	Serum or Plasma (Heparin)	
Fetal Blood (in Maternal Circulation)	Not Available - See Note.	
Fetal Lung Maturation	Amniotic fluid	
Fetoprotein	Amniotic fluid	α-Fetoprotein
Fetoprotein	Fluid	α-Fetoprotein
Fetoprotein	Serum	α-Fetoprotein
Fibrin Degradation Products	Plasma (Citrate)	D-Dimer
Fibrin Split Products	Plasma (Citrate)	D-Dimer
Fibrinogen	Plasma (Citrate)	Fibrinogen: Functional
Fibrinogen	Plasma (Sodium Citrate)	Fibrinogen: Immunological
Fibrinogen: Functional	Plasma (Citrate)	
Fibrinogen: Immunological	Plasma (Sodium Citrate)	
Fiorinal	Serum	Butalbital
Fletcher Factor	Plasma (Citrate)	Prekallikrein
Flexeril	Urine (Random)	Cyclobenzaprine
Fluids	Fluid	Cytology
Fluids	Variable	Pathology Services
Flunitrazepam	Urine (random)	Benzodiazepine Identification
Flunitrazepam	Urine (Random)	
Fluoride	Serum or Plasma (EDTA)	
Fluorocytosine	Not Available - See Note.	5-Fluorocytosine
Fluoxetine	Serum	
Flurazepam	Not Available - See Note.	
Fluvoxamine	Serum	
Folate	Serum	
Folate	Whole blood (EDTA)	
Folic Acid	Serum	Folate
Folic Acid	Whole blood (EDTA)	Folate
Follicle Stimulating Hormone	Serum or Plasma (Hep. or EDTA)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Follitropin	Serum or Plasma (Hep. or EDTA)	Follicle Stimulating Hormone
Fractionated Metanephrines	Plasma (EDTA)	Metanephrines
Fractionated Metanephrines	Urine (24h)	Metanephrines
Fragile X Chromosome	Various	
Free Cortisol	Urine (24h)	Cortisol
Free Digoxin	Serum	Digoxin: Free
Free Glycerol	Not Available - See Note.	Glycerol-Free
Free Hemoglobin	Serum or Plasma (EDTA)	Haptoglobin Electrophoresis
Free Kappa	Serum	Immunoglobulins: Free Light Chains
Free Lambda	Serum	Immunoglobulins: Free Light Chains
Free Light Chains	Urine (Random)	Immunoglobulins: Bence Jones Protein
Free Metanephrines	Plasma (EDTA)	Metanephrines
Free PSA	Serum or Plasma (Li-Hep./EDTA)	Prostate Specific Antigen: Free
Free Protein S	Plasma (Citrate)	Protein S: Free
Free Protoporphyrin	Whole blood (EDTA)	Porphyrins: Quantitation
Free T3	Serum or Plasma (Li-Hep./EDTA)	Triiodothyronine: Free
Free T4	Serum or Plasma (Hep. or EDTA)	Thyroxine: Free
Free Testosterone	Serum	Testosterone: Bioavailable
Free Testosterone	Not Available - See Note.	Testosterone: Free
Free Thyroxin Index	Not Available - See Note.	T3 Resin Uptake
Free Triiodothyronine	Serum or Plasma (Li-Hep./EDTA)	Triiodothyronine: Free
Free Valproic Acid	Serum or Plasma (Heparin)	Valproic Acid: Free
Free beta subunit	Serum or Plasma (Hep. or EDTA)	Chorionic Gonadotrophin
Frisium	Gastric lavage	Clobazam
Frisium	Urine (Random)	Clobazam
Frisium	Serum	Clobazam
Fructosamine	Serum	
Fructose	Semen	
Functional Pre Stress	Plasma (Citrate)	Plasminogen Activator Inhibitor
Furosemide	Not Available - See Note.	
G1PUT	Whole blood (Heparin)	Galactose-1-Phos. Urid. Trans.: Quant
G1PUT	Whole blood (Heparin)	Galactose-1-Phos. Urid. Trans.: Screen
G20210A	Whole Blood (ACD)	Prothrombin Mutation
G6PD	Whole Blood (Heparin)	Glucose-6-Phosphate Dehydrogenase
GBM Antibodies	Serum	Glomerular Base. Membrane



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
		Antibody
GBP-28	Serum	Adiponectin
GGT	Serum	γ-Glutamyltransferase
GH	Serum	Growth Hormone
GHB	Urine (random)	γ-Hydroxybutyrate
GM1 Ganglioside Antibody	Serum	
Gabapentin	Serum or Plasma (Heparin)	
Galactose-1-Phos. Urid. Trans.: Quant	Whole blood (Heparin)	
Galactose-1-Phos. Urid. Trans.: Screen	Whole blood (Heparin)	
Gamma Glutamyl Transferase	Serum	γ-Glutamyltransferase
Gamma hydroxybutyrate	Urine (random)	γ-Hydroxybutyrate
Garamycin	Serum	Gentamicin: Peak
Garamycin	Serum	Gentamicin: Trough
Gastrin	Serum	
Genetic Analysis	Various	Molecular Genetics
Gentamicin: Peak	Serum	
Gentamicin: Trough	Serum	
Gliadin (Deamidated) Antibodies	Serum	
Gliadin Antibodies	Serum	Celiac Profile
Glibenclamide	Not available - See Note.	
Glomerular Base. Membrane Antibody	Serum	
Glucagon	Plasma (EDTA)	
Glucose	Fluid	
Glucose	Plasma (Fluoride)	
Glucose	Urine (24h)	
Glucose-6-Phosphate Dehydrogenase	Whole Blood (Heparin)	
Glutamic Acid Decarboxylase Antibodies	Serum	
Glutamyltransferase	Serum	γ-Glutamyltransferase
Glutathione Peroxidase	Not Available - See Note.	
Glutethimide	Gastric lavage	
Glutethimide	Serum	
Glutethimide	Urine (Random)	
Glyburide	Not available - See Note.	Glibenclamide
Glycated Hemoglobin	Whole blood (EDTA)	Hemoglobin A1c (Glycated)
Glycerol-Free	Not Available - See Note.	
Glycine	CSF	
Glycoprotein Subunit	Serum	α-Glycoprotein Subunit
Glycosylated Hemoglobin	Whole blood (EDTA)	Hemoglobin A1c (Glycated)
Gold	Serum or Plasma (EDTA)	
Gold (Urine)	Not Available - See Note.	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Gonorrhoea	Multiple	
Goop	Urine (random)	γ-Hydroxybutyrate
Goose Antigen	Serum	Avian Precipitins: Goose
Gravol	Serum	Dimenhydrinate
Gravol	Serum	Diphenhydramine
Growth Factor Binding Protein 1	Not Available - See Note.	Insulin-Like Growth Factor B.P. 1
Growth Factor Binding Protein 3	Serum	Insulin-Like Growth Factor B.P. 3
Growth Hormone	Serum	
H Bodies	Whole blood (EDTA)	Heinz Bodies
H.pylori	Multiple	Helicobacter pylori
HBSAg Antibodies	Multiple	Hepatitis B Surface Antigen Antibody
HBc Antibodies	Multiple	Hepatitis B Core Antibody
HBe Antibodies	Multiple	Hepatitis B e Antibody
HBeAg	Multiple	Hepatitis B e Antigen
HBsAg	Multiple	Hepatitis B Surface Antigen
HCG	Serum or Plasma (Hep. or EDTA)	Chorionic Gonadotrophin
HCG	Serum	Chorionic Gonadotrophin: Tumour Marker
HCG Subunit	Serum	α-Glycoprotein Subunit
HDL/LDL Cholesterol	Plasma (EDTA)	Lipoprotein Fract.: Precipitation
HDL/LDL/VLDL Cholesterol	Plasma (EDTA)	Lipoprotein Fract.: Ultracentrif.
HIT	Plasma (Citrate)	Heparin-PF4 Antibodies (HIT)
HLA B27	Whole blood (EDTA)	Human Leukocyte Antigen B27
HLA Typing	Various	Human Leukocyte Antigens Typing
HMB	Whole blood (Heparin)	Porphobilinogen Deaminase
HSV	Multiple	Herpes Simplex Virus
HVA	Urine (24h)	Homovanillate
HVA	Urine (random)	Homovanillate
Hageman Factor	Plasma (Citrate)	Coagulation Factor XII
Halcion	Not Available - See Note.	Triazolam
Haldol	Plasma (Heparin)	Haloperidol
Haloperidol	Plasma (Heparin)	
Haptoglobin	Serum	
Haptoglobin Electrophoresis	Serum or Plasma (EDTA)	
Hb	Plasma (Heparin)	Hemoglobin
Hb Electrophoresis	Whole blood (ACD or EDTA)	Hemoglobin Fractionation
HbA1c	Whole blood (EDTA)	Hemoglobin A1c (Glycated)
Heart Muscle Antibodies	Serum	
Heinz Bodies	Whole blood (EDTA)	
Heinz Body Preparation	Whole blood (EDTA)	Heinz Bodies
Helicobacter pylori	Multiple	
Helper Cells	Whole blood (EDTA)	Lymphocyte Marker Analysis
Hemoglobin	Whole Blood (Heparin)	Methemoglobin
Hemoglobin	Plasma (Heparin)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Hemoglobin A1c (Glycated)	Whole blood (EDTA)	
Hemoglobin Electrophoresis	Whole blood (ACD or EDTA)	Hemoglobin Fractionation
Hemoglobin Fractionation	Whole blood (ACD or EDTA)	
Hemopexin	Serum or Plasma (EDTA)	Haptoglobin Electrophoresis
Hemosiderin	Urine (Random)	
Heparin Co-Factor 2	Plasma (Citrate)	Heparin Co-Factor: Functional
Heparin Co-Factor: Functional	Plasma (Citrate)	
Heparin-Induced Thrombocytopenia Type II	Plasma (Citrate)	Heparin-PF4 Antibodies (HIT)
Heparin-PF4 Antibodies (HIT)	Plasma (Citrate)	
Heparin: Fondaparinux	Plasma (Citrate)	
Heparin: Low Molecular Weight	Plasma (Citrate)	
Heparin: Organ Heparinoid	Plasma (Citrate)	
Heparin: Standard	Plasma (Citrate)	
Hepatic Lipase	Plasma (EDTA)	Post Heparin Lipolytic Activity
Hepatitis A Antibody - IgG	Multiple	
Hepatitis A Antibody - IgM	Multiple	
Hepatitis B Core Antibody	Multiple	
Hepatitis B Surface Antigen	Multiple	
Hepatitis B Surface Antigen Antibody	Multiple	
Hepatitis B e Antibody	Multiple	
Hepatitis B e Antigen	Multiple	
Hepatitis C Antibody	Multiple	
Heptacarboxylic Acid	Faeces (Random)	Porphyryns: Quantitation
Heptacarboxylic Acid	Urine (24h)	Porphyryns: Quantitation
Heptacarboxylic Acid	Urine (Random)	Porphyryns: Quantitation
Heroin: Screen	Urine (Random)	
Herpes Isolation	Multiple	
Herpes Simplex Virus	Multiple	
Heterophile Antibodies	Serum	
Hexacarboxylic Acid	Faeces (Random)	Porphyryns: Quantitation
Hexacarboxylic Acid	Urine (24h)	Porphyryns: Quantitation
Hexacarboxylic Acid	Urine (Random)	Porphyryns: Quantitation
Hexosaminidase	Serum	β -Hexosaminidase
Hg	Urine (24h)	Mercury
Hg	Urine (Random)	Mercury
Hg	Whole blood (EDTA)	Mercury
Hg	Hair	Mercury
Hg	Nail	Mercury
High Molecular Weight Kininogen	Plasma (Citrate)	
High Sensitivity CRP	Serum.	C-Reactive Protein
Histamine	Plasma (EDTA)	
Histamine	Urine (24h)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Histamine	Urine (3h)	
Histone Antibodies	Serum	
Histoplasmosis	Multiple	
Homocysteine	Plasma (EDTA)	
Homogentisate	Urine (Random)	Homogentisic Acid
Homogentisic Acid	Urine (Random)	
Homovanillate	Urine (24h)	
Homovanillate	Urine (random)	
Homovanillic Acid	Urine (24h)	Homovanillate
Homovanillic Acid	Urine (random)	Homovanillate
Human Chorionic Gonadotrophin	Serum or Plasma (Hep. or EDTA)	Chorionic Gonadotrophin
Human Chorionic Gonadotrophin	Serum	Chorionic Gonadotrophin: Tumour Marker
Human Leukocyte Antigen B27	Whole blood (EDTA)	
Human Leukocyte Antigens Typing	Various	
Hydrochlorothiazide	Not Available - See Note.	
Hydrocodone	Urine (Random)	Narcotic Analgesics
Hydromorphone	Urine (Random)	Narcotic Analgesics
Hydroxyamoxapine	Not Available - See Note.	Amoxapine
Hydroxybutyrate	Serum or Plasma (Heparin or Fluoride)	β-Hydroxybutyrate
Hydroxycorticosteroid	Not Available - See Note.	17-Hydroxycorticosteroid
Hydroxymethylbilane Synthase	Whole blood (Heparin)	Porphobilinogen Deaminase
Hydroxyprogesterone	Plasma (Li-heparin)	17-Hydroxyprogesterone
Hydroxyproline (Urine)	Not available - See Note.	
IBC	Serum	Iron
IEP	Not Available - See Note.	Protein Analysis: Immunoelectrophoresis
IFE	Serum	Immunoglobulin D (By IFE)
IFE	Serum	Immunoglobulin E (By IFE)
IFE	Serum	Protein Analysis: Immunofixation
IFE	Not Available - See Note	Protein Analysis: Immunofixation (Urine)
IGD	Serum	Immunoglobulin D
IGF1	Serum	Insulin-Like Growth Factor 1
IGFBP1	Not Available - See Note.	Insulin-Like Growth Factor B.P. 1
IGFBP3	Serum	Insulin-Like Growth Factor B.P. 3
INH	Serum or PL (Hep)	Isoniazid
INH	Urine (Random)	Isoniazid
IRT	Serum	Trypsinogen
Ibuprofen	Serum or plasma (Heparin)	
IgA	Serum	Immunoglobulin A
IgA	Not Available - See Note.	Immunoglobulin A (Saliva)



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
IgA	Not Available - See Note	Immunoglobulin A: CSF
IgA	Serum	Immunoglobulins: A, G and M
IgA	Serum	Protein Analysis: Immunofixation
IgD	Serum	Immunoglobulin D (By IFE)
IgD	Serum	Protein Analysis: Immunofixation
IgE	Serum	Immunoglobulin E
IgE	Serum	Immunoglobulin E (By IFE)
IgE	Serum	Immunoglobulin E Ab. Single Allergen
IgE	Serum	Protein Analysis: Immunofixation
IgG	CSF	Immunoglobulin G
IgG	Serum	Immunoglobulin G
IgG	Serum	Immunoglobulins: A, G and M
IgG	Serum	Protein Analysis: Immunofixation
IgG Precipitins	Serum	Allergic Alveolitis (Lung Serology)
IgG Subclasses	Serum	Immunoglobulin G Subclasses
IgG Subclasses	Serum	Immunoglobulin G1 Subclass
IgG Subclasses	Serum	Immunoglobulin G2 Subclass
IgG Subclasses	Serum	Immunoglobulin G3 Subclass
IgG Subclasses	Serum	Immunoglobulin G4 Subclass
IgG/Albumin Ratio	CSF	
IgG1	Serum	Immunoglobulin G Subclasses
IgG1	Serum	Immunoglobulin G1 Subclass
IgG1	Serum	Immunoglobulin G2 Subclass
IgG1	Serum	Immunoglobulin G3 Subclass
IgG1	Serum	Immunoglobulin G4 Subclass
IgG2	Serum	Immunoglobulin G Subclasses
IgG2	Serum	Immunoglobulin G1 Subclass
IgG2	Serum	Immunoglobulin G2 Subclass
IgG2	Serum	Immunoglobulin G3 Subclass
IgG2	Serum	Immunoglobulin G4 Subclass
IgG3	Serum	Immunoglobulin G Subclasses
IgG3	Serum	Immunoglobulin G1 Subclass
IgG3	Serum	Immunoglobulin G2 Subclass
IgG3	Serum	Immunoglobulin G3 Subclass
IgG3	Serum	Immunoglobulin G4 Subclass
IgG4	Serum	Immunoglobulin G Subclasses
IgG4	Serum	Immunoglobulin G1 Subclass
IgG4	Serum	Immunoglobulin G2 Subclass
IgG4	Serum	Immunoglobulin G3 Subclass
IgG4	Serum	Immunoglobulin G4 Subclass
IgM	CSF	Immunoglobulin M
IgM	Serum	Immunoglobulin M
IgM	Serum	Immunoglobulins: A, G and M
IgM	Serum	Protein Analysis: Immunofixation

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Imipramine	Gastric lavage	
Imipramine	Serum or plasma (EDTA)	
Imipramine	Urine (Random)	
Immune Complexes	Serum	
Immunoglobulin A	Serum	
Immunoglobulin A (Saliva)	Not Available - See Note.	
Immunoglobulin A: CSF	Not Available - See Note	
Immunoglobulin D	Serum	
Immunoglobulin D (By IFE)	Serum	
Immunoglobulin E	Serum	
Immunoglobulin E (By IFE)	Serum	
Immunoglobulin E Ab. Single Allergen	Serum	
Immunoglobulin G	CSF	
Immunoglobulin G	Serum	
Immunoglobulin G Subclasses	Serum	
Immunoglobulin G1 Subclass	Serum	
Immunoglobulin G2 Subclass	Serum	
Immunoglobulin G3 Subclass	Serum	
Immunoglobulin G4 Subclass	Serum	
Immunoglobulin M	CSF	
Immunoglobulin M	Serum	
Immunoglobulins	CSF	IgG/Albumin Ratio
Immunoglobulins: A, G and M	Serum	
Immunoglobulins: Bence Jones Protein	Urine (Random)	
Immunoglobulins: Free Light Chains	Serum	
Immunoreactive Trypsin	Serum	Trypsinogen
Inderal	Not Available - See Note.	Propranolol
Indicans	Urine	
Indoxyl Sulphate	Urine	Indicans
Influenza A	Multiple	
Influenza B	Multiple	
Inhaled Allergens	Serum	
Inhibin A	Serum	
Insulin	Serum or Plasma (Li-Hep./EDTA)	
Insulin Antibodies	Serum	
Insulin-Like Growth Factor 1	Serum	
Insulin-Like Growth Factor B.P. 1	Not Available - See Note.	
Insulin-Like Growth Factor B.P. 3	Serum	
Intralipid	Serum or Plasma (sodium	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
	heparin)	
Intrinsic Factor Antibodies	Serum	
Iodine	Plasma (EDTA)	
Iodine	Urine (24h)	
Iodine	Urine (Random)	
Ionized Calcium	Serum	Calcium: Ionized
Iron	Serum	
Iron	Urine (24h)	
Iron	Tissue	
Islet Cell Antibodies	Serum	Pancreatic Islet Cell Antibodies
Isocoproporphyrin	Not Available - See Note.	
Isoniazid	Serum or PL (Hep)	
Isoniazid	Urine (Random)	
Isonicotinic Acid Hydrazide	Serum or PL (Hep)	Isoniazid
Isopropanol	Gastric lavage	Alcohol Fract. and Quant.
Isopropanol	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Isopropanol	Urine (Random)	Alcohol Fract. and Quant.
Isopropanol	Vitreous Humour	Alcohol Fract. and Quant.
Jo-1	Serum	Extractable Nuclear Antigens Ab.
Jo1	Serum	Extractable Nuclear Antigens Ab.
K	Fluid / Faeces	Potassium
K	Serum or Plasma (Heparin)	Potassium
K	Urine (24h)	Potassium
K	Urine (Random)	Potassium
Kappa Light Chain	Urine (Random)	Immunoglobulins: Bence Jones Protein
Ketalar	Urine (Random)	Ketamine: Screen
Ketamine: Screen	Urine (Random)	
Ketone Bodies	Fluid	Ketones
Ketone Body	Serum or Plasma (Heparin or Fluoride)	β -Hydroxybutyrate
Ketone bodies	Serum or Plasma (Fluoride)	Ketones
Ketones	Fluid	
Ketones	Serum or Plasma (Fluoride)	
Kidney Stone Analysis	Calculus	Calculus Analysis
Killer Cells	Whole blood (EDTA)	Lymphocyte Marker Analysis
L-Lactate	Plasma (Fluoride)	L-Lactic Acid
L-Lactic Acid	Fluid	
L-Lactic Acid	Plasma (Fluoride)	
LAP	Not Available - See note.	Alkaline Phosphatase: Neutrophil
LATS	Serum	Thyroid Receptor Antibody
LCM	Multiple	Lymphocytic Choriomeningitis Virus
LDH	Fluid	Lactate Dehydrogenase
LDH	Serum or plasma (Heparin)	Lactate Dehydrogenase
LH	Serum or Plasma (Hep. or	Luteinizing Hormone



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
	EDTA)	
LKM Antibodies	Serum	Liver Kidney Microsomal Antibodies
LTT	Plasma (Fluoride)	Lactose Absorption Test
La	Serum	Extractable Nuclear Antigens Ab.
Lactate	Serum or plasma (EDTA)	D-Lactic Acid
Lactate	Fluid	L-Lactic Acid
Lactate	Plasma (Fluoride)	L-Lactic Acid
Lactate Dehydrogenase	Fluid	
Lactate Dehydrogenase	Serum or plasma (Heparin)	
Lactose Absorption Test	Plasma (Fluoride)	
Lactose Challenge Test	Plasma (Fluoride)	Lactose Absorption Test
Lactose Tolerance Test	Plasma (Fluoride)	Lactose Absorption Test
Lambda Light Chain	Urine (Random)	Immunoglobulins: Bence Jones Protein
Lamellar Body Count	Amniotic fluid	Fetal Lung Maturation
Lamictal	Serum	Lamotrigine
Lamotrigine	Serum	
Lanoxin	Serum	Digoxin
Largactil	Gastric lavage	Chlorpromazine
Largactil	Serum	Chlorpromazine
Largactil	Urine (Random)	Chlorpromazine
Lead	Fluid	
Lead	Hair	
Lead	Nail	
Lead	Urine (24h)	
Lead	Urine (Random)	
Lead	Whole blood (EDTA)	
Lectopam	Not Available - See Note.	Bromazepam
Legionnaires' Disease	Multiple	
Leptin	Serum or plasma (EDTA)	
Leptospirosis	Multiple	
Leukocyte Alkaline Phosphatase	Not Available - See note.	Alkaline Phosphatase: Neutrophil
Levulose	Semen	Fructose
Librium	Serum	Chlordiazepoxide
Lidocaine	Serum or Plasma (Li Heparin)	
Light Chain	Urine (Random)	Immunoglobulins: Bence Jones Protein
Lipase	Fluid	
Lipase	Serum or plasma (Heparin)	
Lipid Screen	Urine (Random)	Fat: Screen
Lipids	Serum or Plasma (Hep. or EDTA)	Cholesterol
Lipids	Serum or Plasma (Hep.or EDTA)	Triglycerides
Lipolytic Activity	Plasma (EDTA)	Post Heparin Lipolytic Activity



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Lipoprotein Fract.: Precipitation	Plasma (EDTA)	
Lipoprotein Fract.: Ultracentrif.	Plasma (EDTA)	
Lipoprotein Lipase	Plasma (EDTA)	Post Heparin Lipolytic Activity
Lipoprotein Phenotyping	Plasma (EDTA)	Lipoprotein Fract.: Precipitation
Lipoprotein Phenotyping	Plasma (EDTA)	Lipoprotein Fract.: Ultracentrif.
Lipoprotein Ultracentrifugation	Plasma (EDTA)	Lipoprotein Fract.: Ultracentrif.
Lipoprotein a	Serum or plasma (EDTA)	Lipoprotein a
Lipoprotein a	Serum or plasma (EDTA)	
Liquid E	Urine (random)	γ-Hydroxybutyrate
Liquid Ecstasy	Urine (random)	γ-Hydroxybutyrate
Liquid X	Urine (random)	γ-Hydroxybutyrate
Listeriosis	Multiple	
Lithium	Serum	
Liver Cytosol (LC-1) Antibodies	Serum	
Liver Iron	Tissue	Iron
Liver Kidney Microsomal Antibodies	Serum	
Long Acting Thyroid Stimulator	Serum	Thyroid Receptor Antibody
Lorazepam	Urine (random)	Benzodiazepine Identification
Lorazepam	Gastric lavage	
Lorazepam	Serum	
Lorazepam	Urine (Random)	
Loxapac	Not Available - See Note.	Loxapine
Loxapine	Not Available - See Note.	
Lp a	Serum or plasma (EDTA)	Lipoprotein a
Lpa	Serum or plasma (EDTA)	Lipoprotein a
Ludiomil	Gastric lavage	Maprotiline
Ludiomil	Serum or Plasma (EDTA)	Maprotiline
Ludiomil	Urine (Random)	Maprotiline
Luminal	Gastric lavage	Phenobarbital
Luminal	Serum	Phenobarbital
Luminal	Urine (Random)	Phenobarbital
Lupus Anticoagulant	Plasma (Citrate)	Non-Spec. Coagulation Inhib.
Luteinizing Hormone	Serum or Plasma (Hep. or EDTA)	
Lutropin	Serum or Plasma (Hep. or EDTA)	Luteinizing Hormone
Luvox	Serum	Fluvoxamine
Lyme Disease	Multiple	
Lymphocyte Marker Analysis	Whole blood (EDTA)	
Lymphocytic Choriomeningitis Virus	Multiple	
Lysine	Urine (Random)	Cystinuria Monitoring
MDMA	Urine (Random)	Methylenedioxy Amphetamine: Screen



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
MMA	Plasma (EDTA)	Methylmalonic Acid
MMA	Urine (Random)	Methylmalonic Acid
MMF Metabolite	Plasma (EDTA)	Mycophenolic Acid
MPA	Plasma (EDTA)	Mycophenolic Acid
MPO	Serum	Neutrophil Cytoplasmic Antibody
MPO	Serum	Neutrophil Cytoplasmic Antibody (P only)
MPS	Urine (Random)	Mucopolysaccharides: Identification
MPS	Urine (Random)	Mucopolysaccharides: Screen
MS	CSF and 1 mL Serum	Protein Analysis: Oligoclonal Banding
MSS Program	Serum	Maternal Serum Screen
MTHF Reductase	Various	
Macroamylase	Serum	
Macroglobulin	Serum	α-2-Macroglobulin
Macroprolactin	Serum	
Magnesium	Serum	
Magnesium	Urine (24h)	
Magnesium	Urine (Random)	
Magnesium	Whole blood (Heparin)	
Magnesium (Faeces)	Not Available - See Note.	
Manganese	Plasma (EDTA)	
Manganese	Urine (24h)	
Manganese	Urine (Random)	
Manganese	Whole blood (EDTA)	
Maprotiline	Gastric lavage	
Maprotiline	Serum or Plasma (EDTA)	
Maprotiline	Urine (Random)	
Marijuana	Urine (Random)	Cannabinoids: Screen
Maternal Serum Screen	Serum	
Measles	Multiple	
Meat Fibres	Not Available - See Note.	
Mebaral	Not Available - See Note.	Mephobarbital
Mellaril	Gastric lavage	Thioridazine
Mellaril	Serum	Thioridazine
Mellaril	Urine (Random)	Thioridazine
Meperidine	Urine (Random)	Narcotic Analgesics
Meperidine: Screen	Urine (Random)	
Mephenytoin	Not Available - See Note.	
Mephobarbital	Not Available - See Note.	
Meprobamate	Gastric lavage	
Meprobamate	Serum	
Meprobamate	Urine (Random)	
Mercury	Urine (24h)	
Mercury	Urine (Random)	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Mercury	Whole blood (EDTA)	
Mercury	Hair	
Mercury	Nail	
Mesantoin	Not Available - See Note.	Mephenytoin
Mesoporphyrin	Faeces (Random)	Porphyrins: Quantitation
MetHb	Whole Blood (Heparin)	Methemoglobin
Metabolic Screen	Urine (Random)	Cystinuria Monitoring
Metabolic Screen	Urine (Random)	
Metal	Fluid	Aluminium
Metal	Plasma (EDTA)	Aluminium
Metal	Urine (24h)	Aluminium
Metal	Urine (Random)	Aluminium
Metal	Urine (24h)	Antimony
Metal	Urine (Random)	Antimony
Metal	Whole blood (Heparin)	Antimony
Metal	Not Available - See Note.	Antimony - Plasma
Metal	Whole Blood (EDTA)	Arsenic
Metal	Hair	Arsenic
Metal	Nail	Arsenic
Metal	Urine (24h)	Arsenic: Inorganic
Metal	Urine (Random)	Arsenic: Inorganic
Metal	Urine (24h)	Arsenic: Total
Metal	Urine (Random)	Arsenic: Total
Metal	Urine (24h)	Barium
Metal	Urine (Random)	Barium
Metal	Urine (24h)	Beryllium
Metal	Urine (Random)	Beryllium
Metal	Urine (Random)	Bismuth
Metal	Urine (24h)	Boron
Metal	Urine (Random)	Boron
Metal	Erythrocyte	Cadmium
Metal	Urine (24h)	Cadmium
Metal	Urine (Random)	Cadmium
Metal	Whole blood (EDTA)	Cadmium
Metal	Urine (Random)	Calcium
Metal	Erythrocyte	Chromium
Metal	Plasma (EDTA)	Chromium
Metal	Urine (24h)	Chromium
Metal	Urine (Random)	Chromium
Metal	Whole blood (EDTA)	Chromium
Metal	Fluid	Chromium
Metal	Plasma (EDTA)	Cobalt
Metal	Plasma (K2 EDTA)	Cobalt
Metal	Urine (24h)	Cobalt

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Metal	Urine (Random)	Cobalt
Metal	Erythrocyte	Copper
Metal	Fluid	Copper
Metal	Plasma (EDTA)	Copper
Metal	Urine (24h)	Copper
Metal	Urine (Random)	Copper
Metal	Whole Blood (K2 EDTA)	Copper
Metal	Tissue	Copper
Metal	Serum or Plasma (EDTA)	Gold
Metal	Plasma (EDTA)	Iodine
Metal	Urine (24h)	Iodine
Metal	Urine (Random)	Iodine
Metal	Serum	Iron
Metal	Urine (24h)	Iron
Metal	Tissue	Iron
Metal	Fluid	Lead
Metal	Hair	Lead
Metal	Nail	Lead
Metal	Urine (24h)	Lead
Metal	Urine (Random)	Lead
Metal	Whole blood (EDTA)	Lead
Metal	Serum	Magnesium
Metal	Urine (24h)	Magnesium
Metal	Urine (Random)	Magnesium
Metal	Plasma (EDTA)	Manganese
Metal	Urine (24h)	Manganese
Metal	Urine (Random)	Manganese
Metal	Whole blood (EDTA)	Manganese
Metal	Urine (24h)	Mercury
Metal	Urine (Random)	Mercury
Metal	Whole blood (EDTA)	Mercury
Metal	Hair	Mercury
Metal	Nail	Mercury
Metal	Plasma (EDTA)	Molybdenum
Metal	Plasma (K2 EDTA)	Molybdenum
Metal	Urine (24h)	Molybdenum
Metal	Hair	Nickel
Metal	Plasma (EDTA)	Nickel
Metal	Urine (24h)	Nickel
Metal	Urine (Random)	Nickel
Metal	Whole Blood (EDTA)	Nickel
Metal	Fluid / Faeces	Potassium
Metal	Serum or Plasma (Heparin)	Potassium
Metal	Urine (24h)	Potassium

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Metal	Urine (Random)	Potassium
Metal	Plasma (EDTA)	Selenium
Metal	Urine (24h)	Selenium
Metal	Urine (Random)	Selenium
Metal	Erythrocyte	Silver
Metal	Urine (24h)	Silver
Metal	Urine (Random)	Silver
Metal	Urine (24h)	Strontium
Metal	Urine (Random)	Strontium
Metal	Urine (Random)	Sulfur
Metal	Erythrocyte	Thallium
Metal	Urine (24h)	Thallium
Metal	Urine (Random)	Thallium
Metal	Whole Blood (EDTA)	Thallium
Metal	Urine (24h)	Tin
Metal	Urine (Random)	Tin
Metal	Whole Blood (K2 EDTA)	Titanium
Metal	Hair	Uranium
Metal	Urine (24h)	Uranium
Metal	Urine (Random)	Uranium
Metal	Plasma (EDTA)	Vanadium
Metal	Plasma (K2 EDTA)	Vanadium
Metal	Urine (24h)	Vanadium
Metal	Urine (Random)	Vanadium
Metal	Hair	Zinc
Metal	Plasma (EDTA)	Zinc
Metal	Urine (24h)	Zinc
Metal	Urine (Random)	Zinc
Metal	Whole blood (EDTA)	Zinc Protoporphyrin
Metanephrines	Plasma (EDTA)	
Metanephrines	Urine (24h)	
Methadone	Urine (Random)	Narcotic Analgesics
Methadone	Fluid	
Methadone	Serum	
Methadone	Urine (Random)	
Methanol	Gastric lavage	Alcohol Fract. and Quant.
Methanol	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Methanol	Urine (Random)	Alcohol Fract. and Quant.
Methanol	Vitreous Humour	Alcohol Fract. and Quant.
Methaqualone	Urine (Random)	
Methaqualone (Serum)	Not Available - See Note.	
Methemalbumin	Serum or Plasma (EDTA)	Haptoglobin Electrophoresis
Methemalbumin	Plasma (Heparin)	
Methemoglobin	Whole Blood (Heparin)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Methotrexate	Serum	
Methotrimeprazine	Gastric lavage	
Methotrimeprazine	Serum	
Methotrimeprazine: Screen	Urine (Random)	
Methsuximide	Serum or Plasma (Heparin)	
Methylenedioxy Amphetamine: Screen	Urine (Random)	
Methylenetetrahydrofolate Reductase	Various	MTHF Reductase
Methylmalonate	Urine (Random)	Methylmalonic Acid
Methylmalonic Acid	Plasma (EDTA)	
Methylmalonic Acid	Urine (Random)	
Methylphenidate: Screen	Urine (Random)	
Methyprylon	Gastric lavage	
Methyprylon	Serum	
Methyprylon: Urine	Not Available - See Note.	
Mexilitine	Not Available - See Note.	
Mg	Serum	Magnesium
Mg	Urine (24h)	Magnesium
Mg	Urine (Random)	Magnesium
Mg	Whole blood (Heparin)	Magnesium
Microalbumin	Urine (24h)	Albumin
Microalbumin	Urine (Random)	Albumin
Microglobulin	Serum	β -2-Microglobulin
Microglobulin	Urine (Random)	β -2-Microglobulin
micropolyspora	Serum	Allergic Alveolitis (Lung Serology)
Microsomal Antibodies	Serum or Plasma (Na-Heparin)	Thyroid Peroxidase Antibodies
Microsomal Antibody	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Milk Antibodies	Not Available - See Note.	
Mitochondrial Antibodies	Serum	
Mn	Plasma (EDTA)	Manganese
Mn	Urine (24h)	Manganese
Mn	Urine (Random)	Manganese
Mn	Whole blood (EDTA)	Manganese
Mo	Plasma (EDTA)	Molybdenum
Mo	Urine (24h)	Molybdenum
Mogadon	Gastric lavage	Nitrazepam
Mogadon	Serum	Nitrazepam
Mogadon	Urine (Random)	Nitrazepam
Molecular Genetics	Whole Blood (ACD)	Factor V Leiden Mutation
Molecular Genetics	Various	MTHF Reductase
Molecular Genetics	Whole Blood (ACD)	Prothrombin Mutation
Molecular Genetics	Whole Blood (EDTA)	Thiopurine Methyltransferase:



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
		Genotype
Molecular Genetics	Various	
Molybdenum	Plasma (EDTA)	
Molybdenum	Plasma (K2 EDTA)	
Molybdenum	Urine (24h)	
Mono	Serum	Heterophile Antibodies
Monoacetylmorphine	Urine (Random)	Heroin: Screen
Monosialoganglioside GM1 (IgM) Ab.	Serum	GM1 Ganglioside Antibody
Morphine	Urine (Random)	Heroin: Screen
Morphine	Urine (Random)	Narcotic Analgesics
Motrin	Serum or plasma (Heparin)	Ibuprofen
Mucopolysaccharides: Identification	Urine (Random)	
Mucopolysaccharides: Screen	Urine (Random)	
Mullerian Inhibiting Substance	Serum or Plasma (EDTA or Hep)	Anti-Mullerian Hormone
Multimers	Plasma (Citrate)	Von Willebrand's Factor Multimers
Multiple sclerosis	CSF and 1 mL Serum	Protein Analysis: Oligoclonal Banding
Mumps - Soluble	Multiple	
Mumps Virus	Multiple	
Mycophenolic Acid	Plasma (EDTA)	
Mycoplasma Isolation	Not Available - See Note.	
Mycoplasma PN	Multiple	
Myeloperoxidase	Serum	Neutrophil Cytoplasmic Antibody
Myeloperoxidase	Serum	Neutrophil Cytoplasmic Antibody (P only)
Myocardial Antibody	Serum	Heart Muscle Antibodies
Myoglobin	Plasma (Li. Heparin)	
Myoglobin	Urine (Random)	
Mysoline	Gastric lavage	Primidone
Mysoline	Serum	Primidone
N-Acetylprocainamide	Serum	Procainamide
N-Desmethysuximide	Serum or Plasma (Heparin)	Methsuximide
N-Terminal pro Brain Natriuretic Peptide	Serum	B-Type Natriuretic Peptide (NT-Pro)
NDRA	Urine (Random)	Amphetamine
NDRI	Urine (Random)	Methylphenidate: Screen
NH4	Plasma (EDTA)	Ammonia
NMO	Serum	Neuromyelitis Optica Antibody (IgG)
NT-Pro BNP	Serum	B-Type Natriuretic Peptide (NT-Pro)
Na	Fluid / Faeces	Sodium
Na	Urine (24h)	Sodium
Na	Urine (Random)	Sodium
Na	Serum or plasma (Heparin)	Sodium

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Narcotic Analgesics	Urine (Random)	Codeine: Screen
Narcotic Analgesics	Serum	Methadone
Narcotic Analgesics	Urine (Random)	Methadone
Narcotic Analgesics	Urine (Random)	
Narcotic Analgesics	Urine (Random)	Meperidine: Screen
Nebcin	Serum	Tobramycin: Peak
Nebcin	Serum	Tobramycin: Trough
Nembutal	Gastric lavage	Pentobarbital
Nembutal	Serum	Pentobarbital
Nembutal	Urine (Random)	Pentobarbital
Netilmicin	Not Available - See Note.	
Netromycin	Not Available - See Note.	Netilmicin
Neuromyelitis Optica Antibody (IgG)	Serum	
Neuronal Antibodies	Serum	Paraneoplastic Autoantibody Panel
Neuronal Antibody	CSF	Paraneoplastic Autoantibody Panel
Neurontin	Serum or Plasma (Heparin)	Gabapentin
Neutrophil Cytoplasmic Antibody	Serum	
Neutrophil Cytoplasmic Antibody (C only)	Serum	
Neutrophil Cytoplasmic Antibody (P only)	Serum	
Ni	Hair	Nickel
Ni	Plasma (EDTA)	Nickel
Ni	Urine (24h)	Nickel
Ni	Urine (Random)	Nickel
Ni	Whole Blood (EDTA)	Nickel
Nickel	Hair	
Nickel	Plasma (EDTA)	
Nickel	Urine (24h)	
Nickel	Urine (Random)	
Nickel	Whole Blood (EDTA)	
Nicotine: Screen	Urine (Random)	
Nitrazepam	Gastric lavage	
Nitrazepam	Serum	
Nitrazepam	Urine (Random)	
Nitrogen (Feces)	Not Available - See Note	
Nitrogen (Urine)	Not Available - See Note.	
Noludar	Gastric lavage	Methyprylon
Noludar	Serum	Methyprylon
Noludar	Not Available - See Note.	Methyprylon: Urine
Non-Spec. Coagulation Inhib.	Plasma (Citrate)	
Non-Specific Inhibitor	Plasma (Citrate)	Non-Spec. Coagulation Inhib.
Noradrenaline	Plasma (EDTA)	Catecholamines
Noradrenaline	Urine (24h)	Catecholamines



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Norclozapine	Plasma (EDTA or Heparin)	Clozapine
Norepinephrine	Plasma (EDTA)	Catecholamines
Norepinephrine	Urine (24h)	Catecholamines
Norfluoxetine	Serum	Fluoxetine
Normetanephrine	Plasma (EDTA)	Metanephrines
Normetanephrine	Urine (24h)	Metanephrines
Normethsuximide	Serum or Plasma (Heparin)	Methsuximide
Norpace	Plasma (Heparin)	Disopyramide
Norpramin	Gastric lavage	Desipramine
Norpramin	Serum or plasma (EDTA)	Desipramine
Norpramin	Urine (Random)	Desipramine
Nortriptyline	Gastric lavage	
Nortriptyline	Serum or Plasma (EDTA)	
Nortriptyline	Urine (Random)	
Nozinan	Gastric lavage	Methotrimeprazine
Nozinan	Serum	Methotrimeprazine
Nozinan	Urine (Random)	Methotrimeprazine: Screen
Nuclear Antibodies	Fluid	
Nuclear Antibodies	Serum	
Nucleolar Antibody	Serum	Nuclear Antibodies
Nucleotidase	Not Available - See Note.	5' Nucleotidase
Nytol	Serum	Diphenhydramine
OV125	Fluid	CA 125
OV125	Serum	CA 125
Olanzapine	Serum or Plasma (heparin)	
Olanzapine	Urine	
Oligosaccharides	Urine (Random)	
Opiate	Urine (Random)	Codeine: Screen
Opiate	Urine (Random)	Meperidine: Screen
Opiate	Serum	Methadone
Opiate	Urine (Random)	Methadone
Opiate	Urine (Random)	Narcotic Analgesics
Opiates: Screen	Urine (Random)	
Organic Acids	Urine (Random)	
Ornithine	Urine (Random)	Cystinuria Monitoring
Orosomuroid	Serum	α -1-Acid Glycoprotein
Orotic Acid	Urine (Random)	
Osmolality	Fluid / Faeces	
Osmolality	Serum	
Osmolality	Urine (Random or 24h)	
Osmolarity	Fluid / Faeces	Osmolality
Osmolarity	Serum	Osmolality
Osmolarity	Urine (Random or 24h)	Osmolality
Osteocalcin	Serum	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Ovarian Antibodies	Not Available - See Note.	
Oxalate	Urine (24h)	
Oxalate	Urine (Random)	
Oxazepam	Serum	
Oxycodone	Not Available - See Note.	Percocet
Oxycodone	Not Available - See Note.	Percodan
Oxycodone: Screen	Urine (Random)	
P&C ANCA	Serum	Neutrophil Cytoplasmic Antibody
P-ANCA	Serum	Neutrophil Cytoplasmic Antibody (P only)
PAI	Plasma (Citrate)	Plasminogen Activator Inhibitor
PANCA	Serum	Neutrophil Cytoplasmic Antibody
PAP	Variable	Pathology Services
PBG	Urine (Random)	Porphyrin Precursors (ALA & PBG)
PBG Deaminase	Whole blood (Heparin)	Porphobilinogen Deaminase
PCA-1	CSF	Paraneoplastic Autoantibody Panel
PCA-1	Serum	Paraneoplastic Autoantibody Panel
PCP	Urine (Random)	Phencyclidine: Screen
PEP	CSF	Protein Analysis: Electrophoresis
PEP	Serum	Protein Analysis: Electrophoresis
PEP	Urine (Random)	Protein Analysis: Electrophoresis
PEP	CSF and 1 mL Serum	Protein Analysis: Oligoclonal Banding
PHLA	Plasma (EDTA)	Post Heparin Lipolytic Activity
PK	Whole blood (EDTA)	Pyruvate Kinase
PK Screen	Not Available - See Note.	Pyruvate Kinase Screen
PLAP	Serum	Alkaline Phosphatase: Isoenzymes
PO4	Serum	Phosphate
PO4	Urine (24h)	Phosphate
PR3	Serum	Neutrophil Cytoplasmic Antibody
PR3	Serum	Neutrophil Cytoplasmic Antibody (C only)
PSA	Serum	Prostate Specific Antigen
PSA	Serum	Prostate Specific Antigen: Sensitive
PSA-Free	Serum or Plasma (Li-Hep./EDTA)	Prostate Specific Antigen: Free
PSA-Total	Serum	Prostate Specific Antigen
PT	Whole blood (Sodium Citrate)	Prothrombin Time
PTH	Plasma (EDTA) or Serum	Parathyroid Hormone (Intact)
PTH Related Peptide	Serum	Parathyroid Hormone-Related Peptide
PTHrP	Serum	Parathyroid Hormone-Related Peptide
PTT	Whole blood (Citrate)	Partial Thrombo. Time (PTT)
Pancreatic Islet Cell Antibodies	Serum	
Pancreatic lipase	Serum or plasma (Heparin)	Lipase
Parainfluenza 1	Multiple	
Parainfluenza 2	Multiple	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Parainfluenza 3	Multiple	
Paraneoplastic Autoantibody Panel	CSF	
Paraneoplastic Autoantibody Panel	Serum	
Parathyroid Hormone (Intact)	Plasma (EDTA) or Serum	
Parathyroid Hormone-Related Peptide	Serum	
Parathyroid Secretory Protein 1	Plasma (EDTA)	Chromogranin A
Parietal Antibodies	Serum	
Parotid Antibodies	Not Available - See Note.	
Paroxetine	Serum or Plasma (Heparin)	
Parrot Antigen	Serum	Avian Precipitins: Parrot
Partial Thrombo. Time (PTT)	Whole blood (Citrate)	
Pathology Services	Variable	
Paul Bunnell-David	Multiple	
Paxil	Serum or Plasma (Heparin)	Paroxetine
Pb	Fluid	Lead
Pb	Hair	Lead
Pb	Nail	Lead
Pb	Urine (24h)	Lead
Pb	Urine (Random)	Lead
Pb	Whole blood (EDTA)	Lead
Pemphigoid Antibodies	Serum	Skin Antibodies
Pemphigus Antibodies	Serum	Skin Antibodies
Pentacarboxylic Acid	Faeces (Random)	Porphyrins: Quantitation
Pentacarboxylic Acid	Urine (24h)	Porphyrins: Quantitation
Pentacarboxylic Acid	Urine (Random)	Porphyrins: Quantitation
Pentazocine	Urine (Random)	Narcotic Analgesics
Pentobarbital	Gastric lavage	
Pentobarbital	Serum	
Pentobarbital	Urine (Random)	
Pentobarbital	Serum	Thiopental
Pentothal	Serum	Thiopental
Percocet	Urine (Random)	Oxycodone: Screen
Percocet	Not Available - See Note.	
Percodan	Urine (Random)	Oxycodone: Screen
Percodan	Not Available - See Note.	
Perinuclear	Serum	Neutrophil Cytoplasmic Antibody
Perinuclear	Serum	Neutrophil Cytoplasmic Antibody (P only)
Peroxidase Antibodies	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Pertofrane	Gastric lavage	Desipramine
Pertofrane	Serum or plasma (EDTA)	Desipramine



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Pertofrane	Urine (Random)	Desipramine
Pertussis	Multiple	
Pethidine	Urine (Random)	Narcotic Analgesics
Pethidine Hydrochloride	Urine (Random)	Meperidine: Screen
Phadiatop	Serum	Inhaled Allergens
Phenacetin	Not Available - See Note.	
Phencyclidine: Screen	Urine (Random)	
Phenobarbital	Gastric lavage	
Phenobarbital	Serum	
Phenobarbital	Urine (Random)	
Phenolphthalein	Not Available - See Note.	
Phenothiazines: Screen	Urine (Random)	
Phensuximide	Not Available - See Note.	
Phenytoin	Serum	
Phenytoin: Free	Serum	
Phosphate	Serum	
Phosphate	Urine (24h)	
Phospholipid Antibody	Serum	Cardiolipin Antibodies
Phospholipid Antibody	Plasma (Citrate)	Non-Spec. Coagulation Inhib.
Phytanic Acid	Serum or Plasma (Heparin)	
Pigeon	Serum	Allergic Alveolitis (Lung Serology)
Pigeon Antigen	Serum	Avian precipitins: Pigeon
Pinworm Prep	Multiple	
Placidyl	Not Available - See Note.	Ethchlorvynol
Plasma Thromboplastin Antecedent	Plasma (Citrate)	Coagulation Factor XI
Plasmin Inhibitor	Plasma (Citrate)	α -2-Antiplasmin
Plasminogen	Plasma (Citrate)	
Plasminogen Activator Inhibitor	Plasma (Citrate)	
Platelet Antibodies	Not Available - See Note.	Platelet Associated IgG
Platelet Associated IgG	Not Available - See Note.	
Poliovirus 1	Multiple	
Poliovirus 2	Multiple	
Poliovirus 3	Multiple	
Polyoma Virus	Plasma (EDTA)	BK Virus
Polyoma Virus	Urine (Random)	BK Virus
Porphobilinogen	Urine (24h)	Porphyrin Precursors (ALA & PBG)
Porphobilinogen	Urine (Random)	Porphyrin Precursors (ALA & PBG)
Porphobilinogen Deaminase	Whole blood (Heparin)	
Porphyrin Precursors (ALA & PBG)	Urine (24h)	
Porphyrin Precursors (ALA & PBG)	Urine (Random)	
Porphyryns: Quantitation	Faeces (Random)	
Porphyryns: Quantitation	Urine (24h)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Porphyrins: Quantitation	Urine (Random)	
Porphyrins: Quantitation	Whole blood (EDTA)	
Porphyrins: Quantitation - Plasma	Not Available - See Note	
Porphyrins: Screen	Faeces (Random)	
Porphyrins: Screen	Urine (24h)	
Porphyrins: Screen	Urine (Random)	
Post Heparin Lipolytic Activity	Plasma (EDTA)	
Potassium	Fluid / Faeces	
Potassium	Serum or Plasma (Heparin)	
Potassium	Urine (24h)	
Potassium	Urine (Random)	
Powassan Virus	Multiple	
Prealbumin	Serum	
Pregnancy Test	Serum or Plasma (Hep. or EDTA)	Chorionic Gonadotrophin
Prekallikrein	Plasma (Citrate)	
Primidone	Gastric lavage	
Primidone	Serum	
Pro-BNP	Serum	B-Type Natriuretic Peptide (NT-Pro)
ProBNP	Serum	B-Type Natriuretic Peptide (NT-Pro)
Proaccelerin	Plasma (Citrate)	Coagulation Factor V
Procainamide	Serum	
Prochlorperazine	Not Available - See Note.	
Proconvertin	Plasma (Citrate)	Coagulation Factor VII
Procyclidine	Not Available - See Note.	
Progesterone	Plasma (Li-Heparin)	
Prograf	Whole Blood (EDTA)	Tacrolimus (FK506)
Proinsulin	Serum	
Prolactin	Serum	
Proliferating Cell Nuclear Antibodies	Serum	
Pronestyl	Serum	Procainamide
Propafenone	Serum or Plasma (Heparin)	
Propoxyphene	Urine (Random)	Narcotic Analgesics
Propoxyphene	Gastric lavage	
Propoxyphene	Urine (Random)	
Propoxyphene	Serum	
Propranolol	Not Available - See Note.	
Prostate Specific Antigen	Serum	
Prostate Specific Antigen: Free	Serum or Plasma (Li-Hep./EDTA)	
Prostate Specific Antigen: Sensitive	Serum	
Protein	Serum	Immunoglobulin D (By IFE)

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Protein	Serum	Immunoglobulin E (By IFE)
Protein	CSF	
Protein	Fluid	
Protein	Serum	
Protein	Urine (24h)	
Protein	Not Available - See Note.	Protein Analysis: Immunoelectrophoresis
Protein	Serum	Protein Analysis: Immunofixation
Protein	Not Available - See Note	Protein Analysis: Immunofixation (Urine)
Protein Analysis: Electrophoresis	CSF	
Protein Analysis: Electrophoresis	Serum	
Protein Analysis: Electrophoresis	Urine (Random)	
Protein Analysis: Immunoelectrophoresis	Not Available - See Note.	
Protein Analysis: Immunofixation	Serum	
Protein Analysis: Immunofixation (Urine)	Not Available - See Note	
Protein Analysis: Oligoclonal Banding	CSF and 1 mL Serum	
Protein C Activated PTT	Plasma (Citrate)	Activated Protein C Resistance
Protein C Factor Deficiency	Plasma (Citrate)	Activated Protein C Resistance
Protein C: Activity	Plasma (Citrate)	
Protein C: Antigen	Plasma (Citrate)	
Protein C: Functional	Plasma (Citrate)	Protein C: Activity
Protein C: Immunological	Plasma (Citrate)	Protein C: Antigen
Protein Deposition In Tissue	Tissue	
Protein Electrophoresis	Urine (Random)	Immunoglobulins: Bence Jones Protein
Protein Electrophoresis	CSF	Protein Analysis: Electrophoresis
Protein Electrophoresis	Serum	Protein Analysis: Electrophoresis
Protein Electrophoresis	Urine (Random)	Protein Analysis: Electrophoresis
Protein S: Free	Plasma (Citrate)	
Protein S: Free & Total	Plasma (Citrate)	
Prothrombin	Plasma (Citrate)	Coagulation Factor II: Biological
Prothrombin	Plasma (Citrate)	Coagulation Factor II: Ecarin
Prothrombin Mutation	Whole Blood (ACD)	
Prothrombin Time	Whole blood (Sodium Citrate)	
Protoporphyrin	Faeces (Random)	Porphyryns: Quantitation
Protoporphyrin - Free	Whole blood (EDTA)	Porphyryns: Quantitation
Protriptyline	Gastric lavage	
Protriptyline	Serum or plasma (EDTA)	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Protriptyline	Urine (Random)	
Prozac	Serum	Fluoxetine
Pseudocholinesterase	Serum	Cholinesterase: Phenotype
Pseudocholinesterase	Serum	Cholinesterase: Total Activity
Pteroylglutamic Acid	Serum	Folate
Pteroylglutamic Acid	Whole blood (EDTA)	Folate
Purkinje Cell Cytoplasmic Antibody	CSF	Paraneoplastic Autoantibody Panel
Purkinje cell Cytoplasmic Antibodies	Serum	Paraneoplastic Autoantibody Panel
Pyridoxal Phosphate	Plasma (EDTA)	Vitamin B6
Pyridoxine	Plasma (EDTA)	Vitamin B6
Pyruvate	CSF	
Pyruvate	Whole Blood	
Pyruvate Kinase	Whole blood (EDTA)	
Pyruvate Kinase Screen	Not Available - See Note.	
Pyruvic Acid	Whole Blood	Pyruvate
Q Fever	Multiple	
Quaalude	Urine (Random)	Methaqualone
Quaalude	Not Available - See Note.	Methaqualone (Serum)
Quetiapine	Serum	
Quinidine	Serum	
Quinine: Screen	Urine (Random)	
R-Time	Plasma (Citrate)	Reptilase Time
RAST	Serum	Immunoglobulin E Ab. Single Allergen
RBC Cholinesterase	Not Available - See Note.	Acetylcholinesterase: RBC
RBC Folate	Whole blood (EDTA)	Folate
RF	Fluid	Rheumatoid Factor
RF	Serum	Rheumatoid Factor
RMSF	Multiple	Rocky Mountain Spotted Fever
RNP	Serum	Extractable Nuclear Antigens Ab.
Rabies - RFFIT	Multiple	
Rave	Urine (Random)	Methylenedioxy Amphetamine: Screen
Reducing Substances - Faeces	Not available - See Note.	
Reducing Substances - Urine	Not available - See Note.	
Renin	Plasma (EDTA)	
Reptilase Time	Plasma (Citrate)	
Respiratory Syncytial	Multiple	
Restoril	Not Available - See Note.	Temazepam
Reticulin Antibodies	Serum	
Retinol	Serum	Vitamin A
Reverse T3	Not Available - See Note	
Reverse Triiodothyronine	Not Available - See Note	Reverse T3
Rheumatoid Factor	Fluid	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Rheumatoid Factor	Serum	
Rhythmodan	Plasma (Heparin)	Disopyramide
Ribosomal P Antibodies	Serum	
Risperdal	Serum or Plasma (heparin)	Risperidone
Risperidone	Serum or Plasma (heparin)	
Ristocetin Co-Factor Activity	Plasma (Citrate)	Von Willebrand Factor Activity
Ritalin	Urine (Random)	Methylphenidate: Screen
Rivotril	Gastric lavage	Clonazepam
Rivotril	Serum	Clonazepam
Rivotril	Urine (Random)	Clonazepam
Ro	Serum	Extractable Nuclear Antigens Ab.
Rocky Mountain Spotted Fever	Multiple	
Rohypnol	Urine (Random)	Flunitrazepam
rT3	Not Available - See Note	Reverse T3
Rubella	Multiple	
Rubella Screen - Prenatal	Multiple	
Rythmol	Serum or Plasma (Heparin)	Propafenone
S-Sulfocysteine	Urine (Random)	
S.L.E.	Multiple	St Louis Encephalitis Virus
SGOT	Serum	Aspartate Aminotransferase
SGPT	Serum	Alanine Aminotransferase
SHBG	Serum	Sex Hormone Binding Globulin
SNRI	Urine (Random)	Venlafaxine: Screen
SSA	Serum	Extractable Nuclear Antigens Ab.
SSB	Serum	Extractable Nuclear Antigens Ab.
SSRI	Serum	Fluoxetine
SSRI	Serum	Fluvoxamine
SSRI	Serum or Plasma (Heparin)	Paroxetine
SSRI	Serum or plasma (Heparin)	Sertraline
Saccharomyces cer. Antibodies	Serum	
Salicylate	Not Available - See Note.	Percodan
Salicylates	Gastric lavage	Acetylsalicylic Acid
Salicylates	Serum	Acetylsalicylic Acid
Salicylates	Not Available - See Note.	Acetylsalicylic Acid (Urine)
Salivary Cortisol	Saliva	Cortisol
Salivary Duct Antibodies	Not Available - See Note.	Parotid Antibodies
Sandimmune	Whole Blood (EDTA)	Cyclosporine
Sb	Urine (24h)	Antimony
Sb	Urine (Random)	Antimony
Sb	Whole blood (Heparin)	Antimony
Sb	Not Available - See Note.	Antimony - Plasma
Schilling's Test	Not Available - See Note.	
Sci70	Serum	Extractable Nuclear Antigens Ab.
Se	Plasma (EDTA)	Selenium



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Se	Urine (24h)	Selenium
Se	Urine (Random)	Selenium
Secobarbital	Gastric lavage	
Secobarbital	Serum	
Secobarbital	Urine (Random)	
Seconal	Gastric lavage	Secobarbital
Seconal	Serum	Secobarbital
Seconal	Urine (Random)	Secobarbital
Secretory IgA	Not Available - See Note.	Immunoglobulin A (Saliva)
Selenium	Plasma (EDTA)	
Selenium	Urine (24h)	
Selenium	Urine (Random)	
Sensitive Albumin	Urine (24h)	Albumin
Sensitive Albumin	Urine (Random)	Albumin
Sensitive CRP	Serum.	C-Reactive Protein
Sensitive PSA	Serum	Prostate Specific Antigen: Sensitive
Serax	Serum	Oxazepam
Seroquel	Serum	Quetiapine
Serotonin	Serum	
Serotonin	Urine (24h)	
Sertraline	Serum or plasma (Heparin)	
Sex Hormone Binding Globulin	Serum	
Sickle Cell Test	Whole blood (ACD or EDTA)	
Sickling Test	Whole blood (ACD or EDTA)	Sickle Cell Test
Siderophilin	Serum	Transferrin
Silver	Erythrocyte	
Silver	Urine (24h)	
Silver	Urine (Random)	
Silver (Plasma)	Not available - See Note.	
Sinequan	Gastric lavage	Doxepin
Sinequan	Serum or plasma (EDTA)	Doxepin
Sinequan	Urine (Random)	Doxepin
Skeletal Muscle Antibodies	Serum	
Skin Antibodies	Serum	
Sm	Serum	Extractable Nuclear Antigens Ab.
Smooth Muscle Antibodies	Serum	
Sn	Urine (24h)	Tin
Sn	Urine (Random)	Tin
Sodium	Fluid / Faeces	
Sodium	Urine (24h)	
Sodium	Urine (Random)	
Sodium	Serum or plasma (Heparin)	
Solubility Test	Whole blood (ACD or EDTA)	Sickle Cell Test
Somatomedin C	Serum	Insulin-Like Growth Factor 1

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Somatotropin	Serum	Growth Hormone
Spironolactone: Screen	Urine (Random)	
Sputum	Sputum	Cytology
Sputum	Variable	Pathology Services
Sr	Urine (24h)	Strontium
Sr	Urine (Random)	Strontium
St Louis Encephalitis Virus	Multiple	
Stelazine	Serum	Trifluoperazine
Stone Analysis	Calculus	Calculus Analysis
Streptozyme	Not Available - See Note.	Anti-Streptolysin O Screen
Streptozyme	Serum	Anti-Streptolysin O Titre
Striated Antibodies	Serum	Skeletal Muscle Antibodies
Striated Muscle Antibodies	Serum	Skeletal Muscle Antibodies
Strontium	Urine (24h)	
Strontium	Urine (Random)	
Strychnine: Screen	Urine (Random)	
Stuart-Prower Factor	Plasma (Citrate)	Coagulation Factor X
Subunit a	Plasma (Citrate)	Coagulation Factor XIII: Screen
Subunit b	Plasma (Citrate)	Coagulation Factor XIII: Screen
Sugars	Not available - See Note.	Reducing Substances - Faeces
Sugars	Not available - See Note.	Reducing Substances - Urine
Sulfamethoxazole	Not Available - See Note.	
Sulfatides	Urine (Random)	
Sulfhemoglobin	Not Available - See Note.	
Sulfonamides	Not Available - See Note.	Sulfamethoxazole
Sulfonylurea	Not available - See Note.	Glibenclamide
Sulfur	Urine (Random)	
Sulphamethoxazole	Not Available - See Note.	Sulfamethoxazole
Sulphonamides	Not Available - See Note.	Sulfamethoxazole
Surgical Pathology	Variable	Pathology Services
Surmontil	Gastric lavage	Trimipramine
Surmontil	Serum or Plasma (EDTA)	Trimipramine
Surmontil	Urine (Random)	Trimipramine
Synovial Fluid Examination	Synovial fluid	
T&B Cells	Whole blood (EDTA)	Lymphocyte Marker Analysis
T-Cells	Whole blood (EDTA)	Lymphocyte Marker Analysis
T.B.	Multiple	
T3 Resin Uptake	Not Available - See Note.	
T3 Total	Serum	Triiodothyronine: Total
T3U	Not Available - See Note.	T3 Resin Uptake
T7	Not Available - See Note.	T3 Resin Uptake
TBG	Serum	Thyroxine Binding Globulin
TBII	Serum	Thyroid Receptor Antibody
TCO	CSF	Paraneoplastic Autoantibody Panel



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
TCT	Plasma (Citrate)	Thrombin Time
THC	Urine (Random)	Cannabinoids: Screen
TIBC	Serum	Iron
TPA	Not Available - See Note.	Tissue Plasminogen Activator
TPMT	Whole Blood (EDTA)	Thiopurine Methyltransferase: Genotype
TPO	Serum or Plasma (Na-Heparin)	Thyroid Peroxidase Antibodies
TRAB	Serum	Thyroid Receptor Antibody
TRLI	Serum	Trypsinogen
TSH	Serum or Plasma (Hep. or EDTA)	Thyrotrophin Stimulating Hormone
TSH/FSH/LH Subunit	Serum	α -Glycoprotein Subunit
TSI	Serum	Thyroid Receptor Antibody
Tacrolimus (FK506)	Whole Blood (EDTA)	
Tagamet	Urine (Random)	Cimetidine: Screen
Tay-Sach's Disease Marker	Serum	β -Hexosaminidase
TeBG	Serum	Sex Hormone Binding Globulin
Tegretol	Gastric lavage	Carbamazepine
Tegretol	Serum	Carbamazepine
Tegretol Epoxide	Serum	Carbamazepine Epoxide
Telopeptide-N	Urine (Random)	
Temazepam	Not Available - See Note.	
Testosterone Binding Globulin	Serum	Sex Hormone Binding Globulin
Testosterone: Bioavailable	Serum	
Testosterone: Free	Not Available - See Note.	
Testosterone: Total	Serum or Plasma (Heparin)	
Tetanus Anti-Toxin	Multiple	
Tetracyclic Antidepressant	Serum or Plasma (EDTA)	Maprotiline
Tetracyclic Antidepressant	Urine (Random)	Maprotiline
Thallium	Erythrocyte	
Thallium	Urine (24h)	
Thallium	Urine (Random)	
Thallium	Whole Blood (EDTA)	
Theophylline	Serum	
Thermoactinomyces	Serum	Allergic Alveolitis (Lung Serology)
Thiamine	Plasma (EDTA or Heparin)	Vitamin B1
Thiocyanate	Serum	
Thiopental	Serum	
Thiopurine Methyltransferase: Genotype	Whole Blood (EDTA)	
Thioridazine	Gastric lavage	
Thioridazine	Serum	
Thioridazine	Urine (Random)	
Thrombin Clotting Time	Plasma (Citrate)	Thrombin Time

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Thrombin Time	Plasma (Citrate)	
Thyrocalcitonin	Serum	Calcitonin
Thyroglobulin	Serum	
Thyroglobulin Antibodies	Serum or Plasma (Na-Heparin)	
Thyroglobulin Antibody	Serum or Plasma (Na-Heparin)	Thyroid Antibodies
Thyroid Antibodies	Serum or Plasma (Na-Heparin)	
Thyroid Antibodies	Serum or Plasma (Na-Heparin)	Thyroid Peroxidase Antibodies
Thyroid Antibody	Serum or Plasma (Na-Heparin)	Thyroglobulin Antibodies
Thyroid Peroxidase Antibodies	Serum or Plasma (Na-Heparin)	
Thyroid Receptor Antibody	Serum	
Thyroid Stimulating Hormone	Serum or Plasma (Hep. or EDTA)	Thyrotrophin Stimulating Hormone
Thyroid Stimulating Immunoglobulin	Serum	Thyroid Receptor Antibody
Thyrotrophin Stimulating Hormone	Serum or Plasma (Hep. or EDTA)	
Thyrotrophin Receptor Antibody	Serum	Thyroid Receptor Antibody
Thyroxine (T4)	Not Available - See Note.	
Thyroxine Binding Globulin	Serum	
Thyroxine: Free	Serum or Plasma (Hep. or EDTA)	
Tin	Urine (24h)	
Tin	Urine (Random)	
Tissue Plasminogen Activator	Not Available - See Note.	
Tissue Transglutaminase Antibodies	Serum	Celiac Profile
Tissue Transglutaminase Antibodies	Serum	
Tissue Transglutaminase Antibody	Serum	Celiac Profile
Tissue typing	Various	Human Leukocyte Antigens Typing
Titanium	Whole Blood (K2 EDTA)	
TI	Erythrocyte	Thallium
TI	Urine (24h)	Thallium
TI	Urine (Random)	Thallium
TI	Whole Blood (EDTA)	Thallium
Tobramycin: Peak	Serum	
Tobramycin: Trough	Serum	
Tocainide	Not Available - See Note.	

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Tofranil	Gastric lavage	Imipramine
Tofranil	Serum or plasma (EDTA)	Imipramine
Tofranil	Urine (Random)	Imipramine
Tolbutamide	Not Available - See Note.	
Tonocard	Not Available - See Note.	Tocainide
Topamax	Serum	Topiramate
Topiramate	Serum	
Total Bilirubin	Serum	Bilirubin: Total
Total Iron Binding Capacity	Serum	Iron
Total PSA	Serum	Prostate Specific Antigen
Toxoplasmosis	Multiple	
Transcobalamin	Not Available - See Note.	
Transferrin	Serum	
Transplant	Whole Blood (EDTA)	Cyclosporine
Tranxene	Serum	Clorazepate
Trazodone	Serum	
Triazolam	Not Available - See Note.	
Tricyclic Antidepressant	Gastric lavage	Amitriptyline
Tricyclic Antidepressant	Serum or plasma (EDTA)	Amitriptyline
Tricyclic Antidepressant	Urine (Random)	Amitriptyline
Tricyclic Antidepressant	Gastric lavage	Clomipramine
Tricyclic Antidepressant	Serum or plasma (EDTA)	Clomipramine
Tricyclic Antidepressant	Gastric lavage	Desipramine
Tricyclic Antidepressant	Serum or plasma (EDTA)	Desipramine
Tricyclic Antidepressant	Urine (Random)	Desipramine
Tricyclic Antidepressant	Serum or plasma (EDTA)	Doxepin
Tricyclic Antidepressant	Urine (Random)	Doxepin
Tricyclic Antidepressant	Serum or plasma (EDTA)	Imipramine
Tricyclic Antidepressant	Urine (Random)	Imipramine
Tricyclic Antidepressant	Gastric lavage	Nortriptyline
Tricyclic Antidepressant	Serum or Plasma (EDTA)	Nortriptyline
Tricyclic Antidepressant	Urine (Random)	Nortriptyline
Tricyclic Antidepressant	Gastric lavage	Trimipramine
Tricyclic Antidepressant	Serum or Plasma (EDTA)	Trimipramine
Tricyclic Antidepressant	Urine (Random)	Trimipramine
Trifluoperazine	Gastric lavage	
Trifluoperazine	Serum	
Trifluoperazine	Urine (Random)	
Triglycerides	Serum or Plasma (Hep.or EDTA)	
Triiodothyronine - Reverse	Not Available - See Note	Reverse T3
Triiodothyronine: Free	Serum or Plasma (Li-Hep./EDTA)	
Triiodothyronine: Total	Serum	
Trimeprazine	Gastric lavage	



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Trimeprazine	Serum	
Trimeprazine	Urine (Random)	
Trimipramine	Gastric lavage	
Trimipramine	Serum or Plasma (EDTA)	
Trimipramine	Urine (Random)	
Triple Marker Screen	Serum	Maternal Serum Screen
Triptil	Gastric lavage	Protriptyline
Triptil	Serum or plasma (EDTA)	Protriptyline
Triptil	Urine (Random)	Protriptyline
Troponin I	Serum or Plasma (Heparin)	
Troponin T	Serum or Plasma (Li Heparin)	
Trypsin	Not Available - See Note.	
Trypsin-like Immunoreactivity	Serum	Trypsinogen
Trypsinogen	Serum	
Tryptase	Serum	
Tularensis	Multiple	
Turkey Antigen	Serum	Avian Precipitins: Turkey
Tylenol	Gastric lavage.	Acetaminophen
Tylenol	Serum	Acetaminophen
Tylenol	Not Available - See Note.	Percocet
Typhus Epidemic	Multiple	
UNICAP	Serum	Immunoglobulin E Ab. Single Allergen
Ultrasensitive PSA	Serum	Prostate Specific Antigen: Sensitive
Unfractionated Heparin	Plasma (Citrate)	Heparin: Standard
Uranium	Hair	
Uranium	Urine (24h)	
Uranium	Urine (Random)	
Urates	Fluid / Faeces	Uric Acid
Urates	Serum	Uric Acid
Urates	Urine (24h)	Uric Acid
Urea	Fluid / Faeces	
Urea	Serum	
Urea	Urine (24h)	
Uric Acid	Fluid / Faeces	
Uric Acid	Serum	
Uric Acid	Urine (24h)	
Uric Acid	Synovial fluid	Synovial Fluid Examination
Urinary Free Cortisol	Urine (24h)	Cortisol
Uroporphyrin	Faeces (Random)	Porphyryns: Quantitation
Uroporphyrin	Urine (24h)	Porphyryns: Quantitation
Uroporphyrin	Urine (Random)	Porphyryns: Quantitation
Uroporphyrinogen 1 Synthetase	Whole blood (Heparin)	Porphobilinogen Deaminase
V	Urine (24h)	Vanadium
V	Urine (Random)	Vanadium



Index & Cross Reference

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
VDRL	Multiple	
VITD	Serum or Plasma (Heparin)	Vitamin D (1,25-Dihydroxy)
VITD	Serum or plasma (Heparin)	Vitamin D (25-Hydroxy)
VMA	Urine (24h)	Vanillylmandelic Acid
VW Factor Antigen	Plasma (Citrate)	Von Willebrand Assessment
VW Factor Antigen	Plasma (Citrate)	Von Willebrand Factor Antigen
VW Factor Antigen	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
VWF Ristocetin	Plasma (Citrate)	Von Willebrand Assessment
VWF Ristocetin Co-Factor	Plasma (Citrate)	Von Willebrand Factor Antigen
VWF Ristocetin Co-Factor	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Vaccinia-Variola Virus	Multiple	
Valium	Serum	Diazepam
Valproic Acid	Serum	
Valproic Acid: Free	Serum or Plasma (Heparin)	
Vanadium	Plasma (EDTA)	
Vanadium	Plasma (K2 EDTA)	
Vanadium	Urine (24h)	
Vanadium	Urine (Random)	
Vancocin	Serum	Vancomycin: Peak
Vancocin	Serum	Vancomycin: Trough
Vancomycin: Peak	Serum	
Vancomycin: Trough	Serum	
Vanillylmandelic Acid	Urine (24h)	
Varicella-Zoster Virus	Multiple	
Venlafaxine: Screen	Urine (Random)	
Very Long Chain Fatty Acids	Serum or plasma (EDTA)	Fatty Acids, Long Chain
Virus Isolation	Multiple	
Virus Serology	Multiple	
Viscosity	Serum	
Vitamin A	Serum	
Vitamin B1	Plasma (EDTA or Heparin)	
Vitamin B12	Serum	
Vitamin B12 Binding Capacity	Not Available - See Note.	Transcobalamin
Vitamin B6	Plasma (EDTA)	
Vitamin C	Plasma (Heparin)	
Vitamin D	Serum or Plasma (Heparin)	Vitamin D (1,25-Dihydroxy)
Vitamin D	Serum or plasma (Heparin)	Vitamin D (25-Hydroxy)
Vitamin D (1,25-Dihydroxy)	Serum or Plasma (Heparin)	
Vitamin D (25-Hydroxy)	Serum or plasma (Heparin)	
Vitamin E	Serum	
Volatiles	Gastric lavage	Alcohol Fract. and Quant.
Volatiles	Serum or Plasma (Fluoride)	Alcohol Fract. and Quant.
Volatiles	Urine (Random)	Alcohol Fract. and Quant.
Volatiles	Vitreous Humour	Alcohol Fract. and Quant.

**Index & Cross Reference**

Version: 01-Apr-2012

Test Name OR Search Element	Specimen Type	Primary Listing For Test Information
Von Willebrand Assessment	Plasma (Citrate)	
Von Willebrand Assessment	Plasma (Citrate)	Von Willebrand Factor Antigen
Von Willebrand Assessment	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Von Willebrand Factor Activity	Plasma (Citrate)	
Von Willebrand Factor Antigen	Plasma (Citrate)	
Von Willebrand Factor Antigen	Plasma (Citrate)	Von Willebrand Factor Antigen
Von Willebrand Factor Antigen	Plasma (Citrate)	Von Willebrand Ristocetin Co-Factor
Von Willebrand Ristocetin Co-Factor	Plasma (Citrate)	
Von Willebrand's Factor Multimers	Plasma (Citrate)	
W.E.E.	Multiple	Western Equine Encephalitis Virus
Warfarin	Plasma (Heparin)	
Water Channel	Serum	Neuromyelitis Optica Antibody (IgG)
Water Deprivation test	Plasma (EDTA)	Antidiuretic Hormone
Western Equine Encephalitis Virus	Multiple	
XTT	Plasma (Fluoride)	Xylose
Xa Inhibitor	Plasma (Citrate)	Heparin: Fondaparinux
Xa Inhibitor	Plasma (Citrate)	Heparin: Low Molecular Weight
Xa Inhibitor	Plasma (Citrate)	Heparin: Organ Heparinoid
Xa Inhibitor	Plasma (Citrate)	Heparin: Standard
Xanax	Not Available - See Note.	Alprazolam
Xylose	Plasma (Fluoride)	
Xylose Absorption Test	Plasma (Fluoride)	Xylose
Xylose Tolerance Test	Plasma (Fluoride)	Xylose
Yersinia	Multiple	
ZPP	Whole blood (EDTA)	Lead
ZPP	Whole blood (EDTA)	Zinc Protoporphyrin
Zarontin	Serum	Ethosuximide
Zinc	Hair	
Zinc	Plasma (EDTA)	
Zinc	Urine (24h)	
Zinc	Urine (Random)	
Zinc Protoporphyrin	Whole blood (EDTA)	Lead
Zinc Protoporphyrin	Whole blood (EDTA)	
Zn	Hair	Zinc
Zn	Plasma (EDTA)	Zinc
Zn	Urine (24h)	Zinc
Zn	Urine (Random)	Zinc
Zoloft	Serum or plasma (Heparin)	Sertraline
Zyprexa	Serum or Plasma (heparin)	Olanzapine
Zyprexa	Urine	Olanzapine